StabiliT[®] MX
Vertebral Augmentation
System
(For use with StabiliT[®] Bone Cement)

**StabiliT® MX Vertebral Augmentation System
(For use with StabiliT® Bone Cement)**

Instructions for Use (English) (EN)

Important Information – Please Read Before Use

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

INDICATIONS

The StabiliT® MX Vertebral Augmentation System is intended for percutaneous delivery of StabiliT® Bone Cement. The StabiliT® Bone Cement is indicated for the treatment of pathological fractures of the vertebrae using a vertebroplasty or kyphoplasty procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

The VertecoR® MidLine Cement Staging Osteotome, VertecoR® MidLine Osteotome (MLO), VertecoR® StraightLine Cement Staging Osteotome 3.0-11.5 and VertecoR® StraightLine Cement Staging Osteotome 3.0 (SLO) are intended for scraping or coring of bone in the spine.

DESCRIPTION

The StabiliT MX Vertebral Augmentation System is a system for the controlled delivery of StabiliT Bone Cement in the treatment of vertebral compression fractures. The StabiliT MX Vertebral Augmentation System consists of up to eight (8) components (Figure 1).



Figure 1. StabiliT MX Vertebral Augmentation System

- StabiliT Delivery Syringe (#1) - Used for delivering bone cement into the vertebral body.
- StabiliT Delivery Instruments (#2, 3, 4) - Used to provide access to the vertebral body and facilitate delivery of PMMA bone cement in orthopedic procedures including vertebroplasty or kyphoplasty procedures.
 - Master Syringe (#2)
 - Locking Delivery Cannula (LDC) (#3)
 - Cannula (3a)
 - Stylet (3b)
 - StabiliT Introducer with diamond tip stylet (#4) - The StabiliT Introducer is used for percutaneous bone access. The device is packaged with a diamond tip stylet and cannula. The Introducer cannula is 3.6mm in outer diameter. See the chart below for the device length.
 - Cannula (4a)
 - Stylet (4b)
- StabiliT Introducer with bevel tip stylet (#5) in some kit configurations (see product label for details) – The StabiliT Introducer with bevel tip is used for percutaneous bone access. The device is packaged with a bevel tip stylet and cannula. The Introducer cannula is 3.6 mm in outer diameter. See the chart below for the device length.
 - Cannula (5a)
 - Stylet (5b)

Device	Working Length
StabiliT Introducer 3.6-11.5	11 cm
StabiliT Introducer 3.6	13 cm

- StabiliT Bone Cement and Saturate Mixing System (#6) -
 - The StabiliT Bone Cement is indicated for the treatment of pathological fractures of the vertebrae using a vertebroplasty or kyphoplasty procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).
 - The Saturate Mixing System is intended for mixing of StabiliT Bone Cement.
 - Funnel (#6a)
 - Cement Syringe (#6b)
 - Filter (#6c)
 - Locking syringe (#6d)
 - Monomer (#6e)
- VertecoR MidLine Cement Staging Osteotome (MLO) (#7) – The MLO is a sterile single use Osteotome device for scraping or coring of bone to create an area in which cement delivery can be initiated in the treatment of vertebral compression fractures. It is to be used with the StabiliT Introducer. The shaft is 3.0 mm in outer diameter. See the chart below for the working lengths.

Device	Working Length
MidLine Cement Staging Osteotome 3.0-11.5	14.5 cm
MidLine Cement Staging Osteotome 3.0	16.5 cm

- VertecoR StraightLine Cement Staging Osteotome (SLO) (#8) in some kit configurations (see product label for details) – The SLO is a sterile single use osteotome device for scraping or coring of bone to create an area in which cement delivery can be initiated in the treatment of vertebral compression fractures and obtain a biopsy. It is packaged with a blunt tip stylet. It is to be used with the StabiliT Introducer. The shaft is 3.0 mm in outer diameter. See the chart below for working lengths.
 - Cannula (8a)
 - Stylet (8b)

Device	Working Length
VertecoR StraightLine Cement Staging Osteotome 3.0-11.5	13.5 cm
VertecoR StraightLine Cement Staging Osteotome 3.0	15.5 cm

HOW SUPPLIED

- All components are provided sterile. These devices are intended for single use only. DO NOT re-sterilize and/or reuse. DO NOT use if package is open or damaged and notify the manufacturer.

CONTRAINDICATIONS

- The use of this product is contraindicated in patients with coagulation disorders, or with severe pulmonary insufficiency.
- The use of this product is contraindicated in patients with a compromise in the posterior column of the vertebral body or the walls of the pedicles.
- The use of PMMA bone cement is contraindicated in the presence of active or incompletely treated infection at the site where the bone cement is to be applied.
- This product should not be used in patients with sensitivity to any of the components of the PMMA bone cement.

WARNINGS

- Thoroughly read the IFUs for each device including the StabiliT Bone Cement and Saturate Mixing System, (if packaged separately or along with this IFU) before use. The IFU must be followed to perform a procedure using the StabiliT MX Vertebral Augmentation System.
- For safe and effective use, this device should only be used by qualified physicians with training in the clinical procedure in which it is being used. The physician should have specific training, experience, and thorough familiarity with the use and application of this product.
- Always use fluoroscopic guidance with radiographic equipment that provides high quality imaging to avoid patient injury. Use appropriate imaging techniques to confirm correct Working Cannula placement (before and during advancement and after removal), absence of damage to surrounding structures, and appropriate location of delivered bone cement. Imaging, such as venography, can be used to assess the ability of the vertebra to contain the delivered bone cement.
- It is essential to maintain strict sterile technique during the procedure and during all phases of handling this product.
- Precise Working Cannula placement is required for this procedure. Incorrect device placement could result in patient injury.
- The Working Cannula (part of the StabiliT Introducer) is not intended for delivering bone cement. Always use the LDC to deliver bone cement to the vertebral body.

- The Introducer Stylet must be in place inside the Working Cannula during use (e.g., insertion, removal, manipulation).
- Removal of the Working Cannula must be performed by rotation and axial motion. DO NOT bend the cannula sideways, patient injury may occur.
- Dispose of used product per local, state and federal blood borne pathogen controls including biohazard sharps container and disposal procedures.
- DO NOT use if package is opened or damaged. All devices are provided sterile. All devices are sterilized using gamma radiation. These devices are intended for single use only. DO NOT re-sterilizer re-use. Reconditioning, refurbishing, repair, modification, or re-sterilization of the device(s) to enable further use is expressly prohibited, as it may result in patient injury including loss of function and/or infection.
- For devices penetrating bone, DO NOT use if dense bone, including traumatic fractures, is encountered. Device damage resulting in patient injury may occur. Breakage of the device may require intervention or retrieval.
- DO NOT use the StraightLine Osteotome or MidLine Osteotome in fractures due to pararenal or prostatic cancer metastasis of the spine.
- DO NOT use the MidLine Osteotome to scrape or core bone in more than one vertebra.

PRECAUTIONS

- Examine all packaging prior to opening. DO NOT use device if damaged or the sterile packaging is breached. Contact the manufacturer if package is opened or damaged.
- Use the device prior to the Use By Date noted on the device packaging.
- Wear safety glasses or a face shield when delivering the bone cement.
- Ensure that all luer-lock connectors are tightened securely. Improperly secured connections could result in disconnection during injection.
- DO NOT insert the StraightLine Osteotome into the Working Cannula if the Stylet is still attached to the StraightLine Osteotome, as it may result in inadequate coring of the bone.
- Exercise caution in cases involving extensive vertebral destruction and significant vertebral collapse (i.e., the vertebral body is less than 1/3 of its original height). Such cases may lead to a technically difficult procedure.

ADVERSE EVENTS

- Serious adverse events, some with fatal outcome, associated with the use of polymethylmethacrylate (PMMA) include:
 - Myocardial infarction
 - Cardiac arrest
 - Cerebrovascular accident
 - Pulmonary embolism
 - Anaphylaxis
 - Diffusion of the bone cement outside the vertebral body: in the peripheral veins (pulmonary embolism), in the epidural plexus (myelopathy, radiculopathy), in the intervertebral disc
- The most frequent adverse reactions reported with PMMA are:
 - Transitory fall in blood pressure
 - Thrombophlebitis
 - Hemorrhage and hematoma
 - Superficial or deep wound infection
 - Bursitis
 - Short-term cardiac irregularities
 - Heterotopic bone formation
- Other potential adverse events reported for PMMA include:
 - Pyrexia
 - Hematuria
 - Dysuria
 - Bladder fistula
 - Transitory worsening of pain due to heat released during polymerization
 - Nerve entrapment and dysphasia due to extrusion of the bone cement beyond its intended application
 - Adhesions and stricture of the ileum due to heat released during polymerization
- Potential adverse events associated with kyphoplasty or vertebroplasty include:
 - Pneumonia
 - Intercostal neuralgia
 - Collapse of a vertebra adjacent to the one injected, due to an osteoporotic disease
 - Pneumothorax
 - Extravasation of bone cement into soft tissue
 - Fracture of a pedicle
 - Rib fracture in patients with diffuse osteopenia, especially during thoracic kyphoplasty procedures, due to the significant downward force exerted during Working Cannula insertion
 - Compression of the spinal cord with paralysis or loss of feeling
- Adverse events potentially associated with use of the StabiliIT Introducer, MidLine Osteotome or StraightLine Osteotomes include:
 - Nerve injury including puncture of the cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis
 - Pulmonary embolism
 - Hemothorax or pneumothorax
 - Infection, including deep or superficial wound infection
 - Unintended puncture wounds including vascular puncture and dural tear
 - Hemorrhage
 - Hematoma
 - Pain

PREPARATION AND USE

- 1) Check packaging for damage prior to placing contents in sterile field.
- 2) Remove product from package using standard sterile technique.
- 3) Check all components for damage. If the pressure gauge needle on the StabiliIT Delivery Syringe is resting outside the "O" box, do not use.
- 4) Mix bone cement per the StabiliIT Bone Cement and Saturate Mixing System IFU.
- 5) Access the vertebral body using the StabiliIT Introducer.
 - a) Under fluoroscopic guidance direct the StabiliIT Introducer down the pedicle of the vertebral body, while checking A/P lateral images to confirm proper device placement.
 - b) Once the StabiliIT Introducer is positioned in the vertebral body, remove the Stylet with a counterclockwise turn, leaving the Working Cannula in place.

- 6) Create an access channel and a cement staging area to stage bone cement delivery using the SLO and/or the MLO.

WARNING: Use imaging guidance and follow the IFU to avoid patient injury.

- a) Use the SLO to create an access channel and a cement staging area in the bone:
 - i) Remove the Stylet from the SLO and set aside prior to inserting the SLO in the Working Cannula.
 - ii) Remove all instrumentation from the Working Cannula.
 - iii) Using imaging guidance, advance the SLO through the Working Cannula, bring it into contact with bone and verify placement of the tip of the SLO at the intended location. The Introducer/SLO stop limits the SLO shaft to extend 15 mm beyond the distal end of the Working Cannula. Always verify placement under imaging guidance.
 - iv) Using imaging guidance, advance the tip of the SLO to the desired location to carefully scrape or core the bone.

Caution: The location of the Working Cannula in the vertebra should be monitored before and during advancement of the SLO.
 - v) Completely seat the SLO against the Working Cannula to completely create the cement staging space. This may require retraction of the Working Cannula.
 - vi) When scraping or coring is complete, using imaging guidance, remove the SLO from the Working Cannula.
 - b) Use the MLO to create an access channel and a cement staging area in the bone:
 - i) The Articulating Tip is the distal portion of the MLO, the directional arrow points in the direction of the articulating tip bending.
 - ii) The Deployment Handle is the rotating portion of the MLO Handle.
 - (1) Turning the Deployment Handle 1/2 turn clockwise will cause the Articulating Tip to Deploy (fully bend).
 - (2) Turning the Deployment Handle counter-clockwise will cause the Articulating Tip to Un-Deploy (straighten).
 - (3) Do not deploy the Deployment Handle greater than 1/2 turn counter clockwise.
 - iii) Ensure the Articulating Tip is fully extended in the straight position prior to insertion into the Working Cannula.
 - iv) Remove the stylet from the Working Cannula.
 - v) Insert the MLO into the Working Cannula until the first laser mark on the shaft is even with the proximal end of the luer on the Working Cannula. Confirm by image guidance that the distal end of the MLO is at the distal end of the Working Cannula before proceeding. When fully inserted into Working Cannula, the shaft extends approximately 26mm beyond the distal end of the Working Cannula.
 - vi) As the Articulating Tip exits the Working Cannula turn the Deployment Handle to bend the Articulating Tip in the direction of the arrow on the MLO Handle.
 - vii) The MLO can be advanced to the desired position using image guidance.
 - viii) Care should be taken at all times to NEVER strike the arms of the Deployment Handle, especially when rotated from its starting position.
 - ix) Using image guidance (and stabilizing the Working Cannula), the MLO can be carefully withdrawn and advanced multiple times to scrape or core bone until the desired cavity (size and location) is created.
 - x) When the Articulating Tip is substantially deployed, the MLO should not be rotated.
 - (1) The rotation limiting mechanism will slip if the large handle is rotated while the Articulating Tip is substantially articulated.

Caution: The location of the Working Cannula in the vertebra should be monitored before, and during advancement of the MLO through the Working Cannula.
 - xi) When cavity creation is complete, using imaging guidance un-deploy the Articulating Tip and return the Deployment Handle to its starting position.

Caution: Un-Deployment of the Articulating Tip should be done slowly and under imaging guidance while carefully monitoring the position of the tip of the device. Un-Deploy the device by turning the Deployment Handle counter clockwise.
 - xii) Under imaging guidance, remove the MLO from the StabiliIT Introducer's Working Cannula.
 - xiii) The location of the StabiliIT Introducer Working Cannula in the vertebra should be monitored and adjusted if necessary after removal of the MLO.
- 7) Prepare the StabiliIT Delivery Syringe.
 - a) Remove stopcock from StabiliIT Delivery Syringe packaging.
 - b) Connect the StabiliIT Delivery Syringe extension tube to the female connector on the stopcock (opposite from the male connector).
 - c) Ensure the male port on the stopcock is set to the closed position.
 - d) Advance the plunger with enough force to completely remove any air present in the syringe.
 - e) Submerge the stopcock in sterile water (or saline).
 - f) Squeeze the trigger on the StabiliIT Delivery Syringe and pull back the handle to fill the Syringe with fluid. Do so until the entire Syringe is filled.
 - g) Invert the StabiliIT Delivery Syringe, squeeze the trigger, and advance the plunger to remove any remaining air in the syringe and extension tube.
 - h) Squeeze the trigger on the StabiliIT Delivery Syringe and pull back fully to aspirate with sterile water (or saline).

Caution: Inspect the StabiliIT Delivery Syringe tubing and stopcock to ensure that there is no air in the system.
 - 8) Assemble system components
 - a) One minute prior to bone cement delivery (see Table 1) remove the Filter and Funnel Assembly. Clean the Cement Syringe of excess bone cement and completely thread the Cement Syringe into the Locking Delivery Cannula.
 - b) Thread the Master Syringe onto the Cement Syringe.

Caution: Ensure Master Syringe is FULLY threaded onto Cement Syringe before proceeding. Not doing so can cause user injury or device malfunction.
 - c) Attach the StabiliIT Delivery Syringe to the Master Syringe. Securely connect the luer of the Master Syringe to the stopcock at the end of the StabiliIT Delivery Syringe extension tube.

Caution: Do not begin bone cement delivery until saturation and preparation time is complete (see Table 1).

- 9) Cement delivery
- Confirm that the StabiliT Delivery Syringe trigger is released to ensure that the plunger is locked in position.
 - Prime the LDC with bone cement by rotating the StabiliT Delivery Syringe handle in the CLOCKWISE direction. Once bone cement exits LDC tip, stop cement flow by squeezing the trigger on the StabiliT Delivery Syringe. Release trigger to lock the plunger in the withdrawn position. Wipe LDC tip clean.
 - Under fluoroscopic guidance, stabilize the Working Cannula and insert the LDC until the rotating wheel contacts the Working Cannula luer. Rotate the LDC wheel to lock the LDC to the Working Cannula.
 - When prepared to deliver bone cement squeeze the StabiliT Delivery Syringe trigger and push the handle forward until resistance is met and release the trigger. Under fluoroscopic guidance, deliver bone cement by rotating the handle in the CLOCKWISE direction.
 - To stop bone cement delivery, squeeze the trigger on the StabiliT Delivery Syringe. Release trigger to lock the plunger in the withdrawn position. To re-engage, squeeze the trigger and push the handle forward until resistance is met, then release trigger. Continue delivering bone cement by rotating the handle in the CLOCKWISE direction.
- Caution:** To protect the threads of the lock release handle, the gauge must indicate 25 ATM or lower before the quick release mechanism is used to stop flow and relieve pressure.
- Caution:** The quick release mechanism will activate (signaled by a clicking sound) if the operator exceeds the maximum pressure for the StabiliT Delivery Syringe. Once this has occurred, the quick release mechanism may disengage at lower pressures during subsequent attempts to increase pressure.
- Caution:** Following completion of bone cement delivery, remove the LDC from the Working Cannula within 1 minute and immediately insert and lock the Stylet in the Working Cannula. If no additional bone cement delivery is required, remove Introducer (Cannula with Stylet).
- Warning:** Removal of the Working Cannula should only be performed after insertion of the Stylet by rotation and axial motion. DO NOT bend Working Cannula sideways, patient injury may occur.

Table 1: Timing of Various Activities at Different Ambient Temperatures

Activity	Approximate Cumulative Time From Initiation of Saturation (minutes)	
	@ 18-19 °C (65-67 °F)	@ 20-23 °C (68-74 °F)
Bone Cement Saturation and Preparation (See Bone Cement IFU)	0-10 minutes	0-5 minutes
Bone Cement Delivery	10-45 minutes (35 minutes working time)	5-40 minutes (35 minutes working time)


Table 2: Effect of Ambient Temperature on Cement

Temperature °F (°C)	Minimum Setting Time (minutes)
81 (27)	35
73 (23)	51
66 (19)	82

STORAGE & HANDLING

Handle with care. Store in original packaging in a clean, cool, and dry location. Avoid exposure to temperature and humidity extremes.

SYMBOL GLOSSARY

	Caution		Use By
	Consult Instructions For Use		Authorized Representative in the European Community
	Sterilized using Irradiation		Keep away from sunlight
	Sterilized using Ethylene Oxide		Keep away from moisture
	Sterilized using Aseptic Processing Techniques		Flammable
	Lot Number		Single Use Device, DO NOT REUSE
	Catalog Number		Store below 25°C
	No Latex		Manufacturer
	Do not use if package is opened or damaged		Diamond tip
	Short		Bevel Tip
	Long		Cannula Gauge
	Device Length in Centimeters		Introducer Gauge/Locking Delivery Cannula Gauge
	Do not Resterilize	QTY:	Quantity

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd.