

INSTRUCTIONS FOR USE:

STATEMENT OF INTENDED USE: The FLO 30 Hemostasis Valve is intended to maintain a fluid-tight seal around interventional and diagnostic devices during use.

CAUTION: Read manufacturer's instructions for the use of catheters, guide wires, and introducers.

INSTRUCTIONS FOR USE:






Inspect the device prior to use to verify that no damage has occurred during shipping.

1. Connect the valve to a catheter or any device with a female Luer fitting.
2. Check operation by closing the round nut by turning it clockwise to obtain a fluid tight seal.
3. In order to open the seal and insert devices through the valve, turn the round nut in a counterclockwise direction until resistance is felt. Advance the guide wire, catheter or other device an appropriate distance into the vasculature. Close the valve around the shaft of the guide wire, catheter, or other interventional device. This forms a fluid-tight seal, yet does not inhibit movement of the device.

WARNING: It is important that the valve be closed tight enough to prevent fluid leakage, yet not so tight as to restrict function.

REUSE PRECAUTION STATEMENT:

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

	Single Use
	Caution: Consult accompanying document
	Sterilized using Ethylene Oxide
	Non-pyrogenic
	Do not use if package is damaged

Rx Only: CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

CE 0086



www.merit.com



Manufacturer:
 Merit Medical Systems, Inc.
 1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
 1-801-253-1600
 U.S.A Customer Service 1-800-356-3748



Authorized Representative:
 Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
 EC Customer Service +31 43 3588222