ja

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
Do not use if  
package is damaged.  
如包装破损,请勿使用

Proofed by: \_\_\_\_\_ Date: \_\_\_\_\_  
Dimensions checked: \_\_\_\_\_ Copy checked: \_\_\_\_\_

DCN 60711-1  
Vernon Hills, IL  
David Knuth  
02-21-17  
362-23307

en

Intended for use in obtaining core biopsy samples from soft tissue such as breast, kidney, liver, lymph nodes and soft tissue masses. Not intended for use in bone.

**Warning:** This needle is not designed for puncturing the skin.

- A. Obturator Handle
- B. Cannula Handle
- C. Specimen Notch

Excessive force or resistance when the needle is extended may result in a bent specimen notch. A bent specimen notch will interfere with needle function and possibly cause tissue tearing. Excessive lateral force applied against the obturator handle of the needle may interfere with the function of the stop mechanism and may result in deeper tissue penetration than is desired with the potential for tissue or organ damage.

**Instructions for Use**

1. Using sterile technique, prepare the biopsy site.
2. Remove packaging. Remove protective sheath.
3. Test for proper function and alignment by sliding the cannula over the obturator several times. If the parts do not move easily, the unit should be discarded.

**Biopsy Procedure (except Breast Biopsy)**

- 4A. Retract the obturator to cover the specimen notch. Hold the cannula handle and insert the needle assembly so that the specimen notch is within the tissue to be biopsied (a).
- 5A. Retract the cannula to expose the specimen notch. **Note:** Do not move the obturator (b).
- 6A. To cut the tissue, quickly advance the cannula handle (c). If another biopsy of the same tissue is required, leave the cannula positioned at the target biopsy site and follow instructions for multiple biopsy samples.

**Instructions for the Breast Biopsy**

- 4B. Retract the obturator to cover the specimen notch. Hold the cannula handle and insert the needle assembly to the edge of the biopsy site (d).
- 5B. Hold cannula handle firmly. Quickly advance the obturator as far as permitted to position the specimen notch in the tissue to be biopsied (e).
- 6B. To cut the tissue, quickly advance the cannula handle (f). If another biopsy of the same tissue is required, leave the cannula positioned at the target biopsy site and follow instructions for multiple biopsy samples.

**For Multiple Biopsy Samples**

- Hold the cannula handle in a fixed position.
- Lift the obturator shaft slightly to clear the stop mechanism.
- Pull the obturator back to remove it from the cannula.
- Remove the sample from the specimen notch.
- To reinsert the obturator into the cannula, align the obturator shaft with the open side of the cannula handle and the specimen notch facing up. Slide the obturator forward until the obturator shaft locks into the handle.
- Follow one of the sampling procedures described to take additional biopsies.

7. Withdraw the needle assembly with the cannula still advanced over the obturator.
8. Advance the obturator to remove the biopsy specimen from the specimen notch.

Do not use if package is damaged.

**Caution:** For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

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Tru-Cut, Carefusion and the Carefusion logo are trademarks or registered trademarks of CareFusion Corporation, or one of its affiliates.
© 2010, 2011, 2013, 2014, 2017, CareFusion Corporation or one of its affiliates. All rights reserved.