en

**Step 1: Coaxial Introducer Placement**

During the pre-procedure set-up, insert the biopsy device through the coaxial introducer needle to verify compatibility. Adjust sliding stopper on coaxial needle to desired calibration level. Lock in place. Insert coaxial and stylet until sliding stopper rests on skin surface. Image coaxial needle to make sure the tip is on the surface of the lesion. Remove inner stylet.

**Step 2: Sample Notch Size**

Select sample notch size by turning the white coaxial bolt on the Temno Biopsy Needle from a maximum of 19 mm to a minimum of 9 mm as shown by the indicator window. The amount of sampling notch exposure and needle advancement are measured from the distal end of the coaxial introducer.

**Step 3: Temno Needle Spring Charge**

Charge Temno needle by pulling back hard on the plunger until a firm click is felt indicating that the needle spring is charged.

**Step 4: Temno Needle Placement**

With the needle stylet fully retracted, insert the Temno needle into the coaxial introducer and gently connect the coaxial's hub onto the Temno's white coaxial bolt until snug.

**Step 5: Capture Biopsy Specimen**

Advance the stylet forward, with thumb on plunger, to deliver the specimen sampling notch into the area to be biopsied. Release the cutting cannula by fully depressing the plunger with thumb. This will cut a tissue core.

**Step 6: Retrieve Biopsy Specimen**

Withdraw Temno needle from coaxial introducer. Pull back on plunger as in Step 3 above until firm click is felt. Carefully advance the stylet forward to expose specimen notch and biopsy core. Remove tissue sample.

**Step 7. Obtaining Multiple Samples**

Additional tissue samples may be obtained by repeating Step 4 through Step 6 above. Note that the distal tip of the coaxial introducer remains proximal to target area.

**Caution:** Test firing the needle may cause damage to the cutting cannula and is not recommended. If test firing is performed, the cutting cannula may be protected by holding the cannula with moderate pressure between the thumb and fingers during firing.

**Note:** The needle stylet will be completely covered by the cutting cannula after charging the device.

**Caution:** For single use only. Re-use may result in a non-functional product or contribute to cross contamination. Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

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Temno is a trademark or registered trademark of CareFusion Corporation, or one of its subsidiaries.
© 2010, 2011, 2014, CareFusion Corporation or one of its subsidiaries. All rights reserved.
**Made in Dominican Republic by CareFusion DR 203 LTD**
Sponsored by **CareFusion Australia 316 Pty. Ltd.**
**CareFusion New Zealand 313 Ltd.**