

# THE SAVI® BREAST BRACHYTHERAPY APPLICATOR: Partial Breast Radiation Therapy for Early-Stage Breast Cancer



## WHAT IS SAVI?

The Strut Adjusted Volume Implant (SAVI) applicator delivers a form of accelerated partial breast radiation (APBI) known as breast brachytherapy, which delivers radiation from inside the breast. This allows physicians to precisely deliver treatment to the tumor cavity and surrounding tissue, and offers two key advantages: it spares healthy tissue from unnecessary radiation and reduces treatment time from 6 weeks to just 2–5 days\*.

## HOW DOES IT WORK?


SAVI uses a bundle of catheters to deliver radiation. The catheters are expanded to fit the size and shape of the cavity, and a tiny radioactive seed travels through each catheter. This allows every catheter to deliver an individual dose of radiation, so physicians can customize treatment for every patient, regardless of the size, shape, or location of their tumor cavity.

## BENEFITS OF SAVI BREAST BRACHYTHERAPY

- Customizes radiation to individual patient
- Excellent local control of cancer
- Reduces toxicity from unnecessary radiation
- Treatment lasts just 2–5 days\*


\*based on clinicians' preferred treatment protocol. Khan, Atif J. et al. Three-Fraction Accelerated Partial Breast Irradiation (Apbi) Delivered With Brachytherapy Applicators Is Feasible And Safe: First Results From The Triumph-T Trial. International Journal of Radiation Oncology. 2019. Online publication.

**1 PLACEMENT**  
The Physician gently **inserts SAVI in a closed position through a small incision**, and then **expands the catheters to conform** to the shape of the cavity. The ends of the catheters remain





outside the breast for radiation delivery. **SAVI remains in the breast** during the entire course of treatment.

**2 CUSTOMIZED TREATMENT PLANNING**  
The radiation oncologist performs a CT scan and uses those images to create a **customized treatment plan to deliver radiation to the areas that**



**need it most**, while protecting the skin, chest wall, heart and lungs.

**3 RADIATION DELIVERY**  
Treatment is **delivered twice a day for 2–5 days**: Each treatment lasts approximately 10 minutes. The ends of the catheters are connected to a device—called an **HDR afterloader**—and a computerized system painlessly delivers a wire with a **tiny radioactive seed into each catheter**.



The seeds are completely removed after each treatment—**no radiation remains in the patient's body** between treatments. After the last treatment, the physician closes the device and removes it through the same incision in which it was inserted.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



merit.com  
ciannamedical.com

©2019 Merit Medical Systems, Inc. All rights reserved.

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
Main: +1.801.253.1600  
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733