

Dr. Deepankar Sharma Director of Interventional Pulmonology Columbus Regional Health, Columbus, IN




PRESENTATION

A 48-year-old male with a medical history significant for diabetes, severe COPD, and 40 pack years of smoking, presented with a progressive cough and shortness of breath over a 3-month duration.

WORK UP

CT chest revealed right lung mass, mediastinal lymphadenopathy and severe stenosis of the bronchus intermedius and the right lower lobe. Imaging also revealed spinal lesions concerning for metastatic lung cancer.



INTERVENTION


Endobronchial mucosal invasion was found in distal bronchus intermedius and right lower lobe on bronhoscopic airway exam. There was moderate stenosis in the bronchus intermedius and severe stenosis in right lower lobe bronchus. Nd:YAG laser was used to treat the endobronchial mucosal tumor. Most of the disease was infiltrative rather than exophytic so no mechanical excision was performed.

The decision was made to place AERO[®] Tracheobronchial stents to preserve the patency of the airway while patient started treatment and showed significant response as his respiratory symptoms were getting worse.


A 12mm x 30mm AERO[®] Tracheobronchial stent was placed in the bronchus intermedius and a 10mm x 10mm AEROmini[®] Tracheobronchial stent in the right lower lobe in a telescoping fashion.

An EBUS TBNA biopsy sample was obtaned from the LN 11L and Station 7, with an on-site cytology preliminary identification of small cell carcinoma.

Bronchus Intermedius and Right Lower Lobe Stenosis




Bronchus Intermedius and Right Lower Lobe Telescoping Stents



FOLLOW UP


The final staging of the small cell lung cancer was determined to be T4N34M1. The patient recieved radiation to L3 lesion and completed Carboplatin, Etoposide and Tecentriq[®]. The patient is being closely followed and is doing well at post 6-months treatment. The future treatment plan included another bronchoscopy to assess for potential distal stent removal, but the stents have not been removed as the patient is feeling great and chose to delay the follow-up bronchoscopy to have the stents removed.

CT-Chest at 6-Month Interval Comparing Right Lower Lobe Airway Stenosis



CONCLUSIONS

The use of AERO® Tracheobronchial stents in a telescoping configuration allows the stenting of a complete airway for treatment and relief of stenosis.



The above information, which may include both Merit and third-party products, are for the practitioner's convenience and for general information purposes only. This information does not constitute medical or legal advice, nor is it meant to endorse or guarantee the suitability of any of the referenced products or methods for any specific patient or procedure. Before using any product, refer to the Instructions for Use (IFU) for indications, contraindications, warnings, precautions, and directions for use.