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- 2006 SIR Guidelines for the Use of Retrievable and Convertible Vena Cava Filters;⁸⁴
- 2003 SIR Clinical Practice Guidelines;⁸⁵
- 2008 European Society of Cardiology (ESC) Guidelines on the Diagnosis and Management of Acute Pulmonary Embolism;⁷⁹
- 2018 ACR Appropriateness Criteria for the Management of Venous Thromboembolism and Inferior Vena Cava Filters;⁸⁶
- 2018 Kidney Disease Outcome Quality Initiative (KDOQI) Clinical Practice Guidelines for Vascular Access.⁸⁷

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Naslov standarda	Skladnost družbe Merit, datum/različica	Skladnost družbe Merit z zadnjim stanjem tehničnega razvoja (popolno/delno/ne, glede na razdelke veljavnih standardov)	Utemeljitev za ne (skladnost) ali delno skladnost
Označevanje			
Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements (ISO 15223-1:2016, Corrected version 2017-03)	EN ISO 15223-1:2016	popolno	N/A
Terminology, Symbols and Information Provided with Medical Devices; Information Supplied by the Manufacturer with Medical Devices	EN 1041:2008+A1:2013	popolno	N/A
Splošni standardi – sterilizacija			
Sterilization of Medical Devices – Validation and Routine Control of Ethylene Oxide Sterilization	EN ISO 11135:2014 & EN ISO 11135:2014/A1:2019	popolno	N/A
Product Adoption and process equivalency for ethylene oxide sterilization	AAMI TIR28:2016	popolno	N/A
Classification of Air Cleanliness, Clean rooms & Associated Controlled Environments. Part 1: Classification of air cleanliness	EN ISO 14644-1:2015	popolno	N/A
Cleanrooms and associated controlled environments – Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration	EN ISO14644-2:2015	popolno	N/A

Naslov standarda	Skladnost družbe Merit, datum/različica	Skladnost družbe Merit z zadnjim stanjem tehničnega razvoja (popolno/delno/ne, glede na razdelke veljavnih standardov)	Utemeljitev za ne (skladnost) ali delno skladnost
Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products	EN ISO 11737-1:2018	popolno	N/A
Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide sterilization residuals	EN ISO 10993-7:2008 & EN ISO 10993-7:2008/AC:2009	popolno	N/A
Bacterial Endotoxins Test	ANSI/AAMI ST72:2019	popolno	N/A
Splošni standardi – sistemi kakovosti			
Quality Systems – Medical Devices – Quality Management Systems. Requirements for Regulatory Purposes	EN ISO 13485:2016	popolno	N/A
Obvladovanje tveganja			
Medical Devices – Application of Risk Management to Medical Devices	EN ISO 14971:2019	popolno	N/A
Biološka varnost			
Biological Evaluation of Medical Devices – Part 1: Evaluation and testing	ISO 10993-1:2018	popolno	N/A
Biological evaluation of medical devices – Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	ISO 10993-3:2014	popolno	N/A
Biological evaluation of medical devices – Part 4: Selection of tests for interactions with blood	ISO 10993-4:2017	popolno	N/A
Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	ISO 10993-5:2009	popolno	N/A
Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization	ISO 10993-10: 2010	popolno	N/A
Biological evaluation of medical devices – Part 11: Tests for systemic toxicity	ISO 10993-11:2017	popolno	N/A
Biological Evaluation of Medical Devices – Part 12: Sample preparation and reference materials	ISO 10993-12:2012	popolno	N/A
Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterization of materials	ISO 10993-18:2020	popolno	N/A
Standard Guide for Biocompatibility of Medical Device Packaging Materials	ISO 10993-19:2020	popolno	N/A

Naslov standarda	Skladnost družbe Merit, datum/različica	Skladnost družbe Merit z zadnjim stanjem tehničnega razvoja (popolno/delno/ne, glede na razdelke veljavnih standardov)	Utemeljitev za ne (skladnost) ali delno skladnost
Klinična ocena			
Clinical evaluation: Guide for manufacturers and notified bodies	MEDDEV 2.7/1 Rev4	popolno	N/A
Nadzor zasnov – kateter			
Intravascular catheters – Sterile and single-use catheters – Part 1: General requirements	ISO 10555-1:2013 & ISO 10555-1:2013/ Amd 1:2017	popolno	Razdelki brez preskušanja: 4.9 Flowrate; 4.10 Power injection; 4.11 Side holes. Za zankasti kateter za vse velja N/A.
Small-bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications	ISO 80369-7:2016	popolno	N/A
Nadzor zasnov – zankasta elektroda			
Sterile Single-Use Intravascular Catheter Introducers	ISO 11070:2014/A1:2018	popolno	V poštev prideta le razdelek 4 “General Requirements” in razdelek 8 “Additional requirements for guidewires”.
Uporabnost			
Medical Devices – Application of usability engineering to medical devices	IEC 62366-1:2015 & IEC 62366-1:2015/ COR1:2016	delno	V skladu z ISO62366-1: 2015 Annex C. Izdelek sproščen za proizvodnjo pred letom 2015, zato velja le ISO62366-1:2015 Annex C.
Embalaza			
Packaging for Terminally Sterilized Medical Devices. Part 1: Requirements for materials, sterile barrier systems, and packaging systems.	EN ISO 11607-1:2020	popolno	N/A
Packaging for Terminally Sterilized Medical Devices. Part 2: Validation requirements for forming, sealing and assembly processes	EN ISO 11607-2:2020	popolno	N/A

Naslov standarda	Skladnost družbe Merit, datum/različica	Skladnost družbe Merit z zadnjim stanjem tehničnega razvoja (popolno/delno/ne, glede na razdelke veljavnih standardov)	Utemeljitev za ne (skladnost) ali delno skladnost
Packaging – Complete, filled transport packages and unit loads – Conditioning for testing	EN ISO 2233:2001	popolno	N/A
Standard Practice for Performance Testing of Shipping Containers and Systems – ASTM D4169 - 16	ASTM D4169 - 16	popolno	N/A
Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurization (Bubble Test) – ASTM F 2096 - 11 (2019)	ASTM F2096 - 11 (2019)	popolno	N/A
Standard Test Method for Detecting Seal Leaks in porous Medical Packaging by Dye Penetration – ASTM F 1929	ASTM F1929 - 15	popolno	N/A
Standard Test Method for Seal Strength of Flexible Barrier Materials – ASTM F88	ASTM F88 / F88M - 15	popolno	N/A
Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices – ASTM F1980	ASTM F1980 - 16	popolno	N/A
Klinično spremljanje po dajanju na trg			
Post Market Clinical Follow-up studies	MEDDEV 2.12/2 Rev2	popolno	N/A
Vigilanca			
Guidelines on a Medical Devices Vigilance System	MEDDEV 2.12/1 Rev8	popolno	N/A

9.0 Reference

1. Alkhouli M, Aljohani S, Kucera R, Moiduddin N, Hull R. A Double-Snare Technique for Safe Retrieval of Embolized Left Atrial Appendage Occluders. *JACC Cardiovasc Interv.* 2017;10(5):e55-e56.
2. Gabriels J, Beldner S, Khan M, Zeitlin J, Jadonath R, Patel A. Embolization of Watchman Device Following a Hybrid Radiofrequency Ablation of Atrial Fibrillation and Watchman Implantation Procedure. *Journal of Cardiovascular Electrophysiology.* 2017;28(7):835-836.
3. Katta N, Gautam S, Webel R. Successful Percutaneous Retrieval of Embolized Septal Occluder Device from Aortic Arch and Placement of a Newer Septal Occluder Device in Combined Procedure. *Case Rep Cardiol.* 2016;2016:1032801.
4. Poomipanit P, Levi D, Shenoda M, Tobis J. Percutaneous retrieval of the locked helex septal occluder. *Catheter Cardiovasc Interv.* 2011;77(6):892-900.
5. Le J, Grigorian A, Chen S, Kuo IJ, Fujitani RM, Kabutey NK. Novel endovascular technique for removal of adherent PICC. *J Vasc Access.* 2016;17(6):e153-e155.
6. Robbins AB, Patel H. Transvascular retrieval of a catheter remnant from the peripheral vein of a preterm neonate. *J Vasc Access.* 2008;9:299-300.

7. Yu VY, Khaja MS, Liles A, Downing T, Williams DM. Endovascular Removal of Fragmented and Embedded Central Venous Catheters via Endoluminal Balloon Traction Method. *J Vasc Interv Radiol*. 2019;30(12):2045-2047 e2041.
8. Anwar S, Little C, Kingston R, Diakar P, Kaikini R. Line stripping: An underused, safe and effective option for blocked tunnelled hemodialysis catheters. *Hemodial Int*. 2013;17(1):141.
9. Beigi AA, Yaribakht M, Sehhat S. Four cases of adhered permanent double lumen hemodialysis catheters (Permcath). *Arch Iran Med*. 2013;16(4):243-245.
10. Desai KR, Laws JL, Salem R, et al. Defining Prolonged Dwell Time: When Are Advanced Inferior Vena Cava Filter Retrieval Techniques Necessary? *Circulation: Cardiovascular Interventions*. 2017;10(6).
11. Desai KR, Xiao N, Karp J, et al. Single-session inferior vena cava filter removal, recanalization, and endovenous reconstruction for chronic ilio caval thrombosis. *J Vasc Surg Venous Lymphat Disord*. 2019;7(2):176-183.
12. Etkin Y, Glaser JD, Nation DA, et al. Retrievable inferior vena cava filters can always be removed using "fall-back" techniques. *J Vasc Surg Venous Lymphat Disord*. 2015;3(4):364-369.
13. Kuo WT, Robertson SW. Bard Denali Inferior Vena Cava Filter Fracture and Embolization Resulting in Cardiac Tamponade: A Device Failure Analysis. *Journal of Vascular and Interventional Radiology*. 2015;26(1):111-115.e111.
14. Lindquist J, Reghunathan A, Brown M, et al. 04:21 PM Abstract No. 198 Comparison of the Captus device versus conventional snare for IVC filter retrieval. *Journal of Vascular and Interventional Radiology*. 2019;30(3):S89-S90.
15. Najafi A, Koulia K, Aubert P, Binkert CA. Multi loop snare technique for difficult inferior vena cava filter retrievals. *CVIR Endovasc*. 2018;1(1):33.
16. Nissim L, Romano W. Successful retrieval from right ventricular wall of an embolized Bard Meridian filter fragment. *J Vasc Interv Radiol*. 2013;24(12):1933-1934.
17. Roh S, Laroia S. Complications during Conversion of VenaTech Inferior Vena Cava Filter. *J Vasc Interv Radiol*. 2017;28(2):304-305.
18. Sinha SK, Mahotra A, Abhishekh NK, et al. Acute Stent Loss and Its Retrieval of a Long, Tapering Morph Stent in a Tortuous, Calcified Lesion. *Cardiol Res*. 2018;9(1):63-67.
19. Ostapenko A, Richard MN, Dietzek AM. Rare Case of Acute Type B Dissection Causing Complete Collapse of EVAR Stent. *Vasc Endovascular Surg*. 2019;53(5):420-423.
20. Havelka GE, Tomita TM, Malaisrie SC, Ho KJ, Eskandari MK. Two Cases of Aortic Intimal Intussusception During Endovascular Repair of an Acute Type B Dissection. *J Endovasc Ther*. 2016;23(3):521-528.
21. DerDerian T, Ascher E, Hingorani A, Jimenez R. A rare complication of a retained wire during endovascular abdominal aortic aneurysm repair. *Ann Vasc Surg*. 2013;27(8):1183 e1111-1185.
22. Moshayedi P, Jadhav AP. Direct Aspiration Catheter Fracture and Retrieval during Neurothrombectomy. *Interv Neurol*. 2018;7(3-4):148-152.
23. Yildiz BS, Alihanoglu YI, Alur I, Evrengul H, Kaya D. A rare complication: an attempt of retrieval of an aortic valve wrapped with pig tail catheter during transcatheter aortic valve implantation. *Cardiovasc Revasc Med*. 2015;16(6):376-377.
24. Kawamura A, Nishiyama N, Kawakami T, Fukuda K. Retrieval of embolized Amplatzer septal occluder using ablation catheter and triple-loop snare. *Cardiovasc Interv Ther*. 2014;29(4):350-353.

25. Funaki B. Catheter Pinch off with Foreign Body Retrieval. *Semin Intervent Radiol.* 2017;34(1):81-84.
26. Huang P-C, Lui D-W, Ng K-K, Ng S-H, Wan Y-L. Percutaneous transvascular retrieval of the dislodged Port-A catheter using EN Snare system via femoral vein approach. *Chin J Radiol.* 2009;34:145-152.
27. Levin R, Valverde IV, Porcile R. Percutaneous removal of knotted Swan-Ganz catheter. *Revista Argentina de Cardiologia.* 2019;87(4):309.
28. Naeim HA, Abuelatta R, Alamodi O, et al. Transcatheter Closure of Aorta to Left Atrium Fistula during Active Prosthetic Valve Endocarditis as a Bridge for Surgery: A Case Report. *CASE (Phila).* 2019;3(6):255-258.
29. Nakabayashi K, Nomura H, Isomura D, Sugiura R, Oka T. Successful Retrieval of a Dismembered Central Venous Catheter Stuck to the Right Pulmonary Artery Using a Stepwise Approach. *Case Reports in Cardiology.* 2016;2016:1-6.
30. Norwa J, Darocha S, Dobosiewicz A, et al. Percutaneous retrieval of a fractured portacath fragment in two patients undergoing long-term chemotherapy. *OnCOReview.* 2016;6(2):A57-61.
31. Smith JC, Black MD, Smith DC. Occlusion balloon catheter-assisted removal of a wedged catheter fragment without a free end from a pulmonary artery. *Cardiovasc Intervent Radiol.* 2009;32(6):1308-1311.
32. Teragawa H, Sueda T, Fujii Y, et al. Endovascular technique using a snare and suture for retrieving a migrated peripherally inserted central catheter in the left pulmonary artery. *World J Cardiol.* 2013;5(9):369-372.
33. Ahmed O, Masrani A, Turba UC, March R, Ginsburg M, Arslan B. Percutaneous Endovascular Retrieval of a Migrated Permanent VenaTech LP Inferior Vena Cava Filter from the Right Ventricle. *Journal of Vascular and Interventional Radiology.* 2017;28(9):1304-1306.
34. Hannawa KK, Good ED, Haft JW, Williams DM. Percutaneous Extraction of Embolized Intracardiac Inferior Vena Cava Filter Struts Using Fused Intracardiac Ultrasound and Electroanatomic Mapping. *J Vasc Interv Radiol.* 2015;26(9):1368-1374.
35. Knavel EM, Woods MA, Kleedehn MG, Ozkan OS, Laeseke PF. Complex Inferior Vena Cava Filter Retrieval Complicated by Migration of Filter Fragment into the Aorta and Subsequent Distal Embolization. *Journal of Vascular and Interventional Radiology.* 2016;27(12):1865-1868.
36. Sharma AK, Ganatra S, Hansen J, et al. A Dual-Snare Percutaneous Retrieval of Venous Stent Embolization to the Right Heart. *JACC Cardiovasc Interv.* 2017;10(12):e111-e113.
37. Sugahara T, Nakashima K, Suzuki R, et al. Percutaneous removal of a fully expanded SMART stent from the pulmonary artery using various adjunctive techniques. *Jpn J Radiol.* 2012;30(8):688-692.
38. Hsu JT, Nakamura S, Yamamoto R, et al. Acute-on-chronic massive pulmonary embolism treated with 3-loop snare. *Circ J.* 2008;72:1909-1913.
39. Sherev D, Azizi K, Azimi NA, Van Nordheim S, Moreno-Cabral R. Delayed Management of Partial Aortic Valve Avulsion After Transcatheter Aortic Valve Replacement. *Ann Thorac Surg.* 2017;104(1):e17-e18.
40. Bowers RW, Beckett D, Paisey JR, Balasubramaniam RN, Sopher SM. Successful release of an entrapped circular mapping catheter using a snare and a multidisciplinary approach. *Kardiol Pol.* 2015;73(6):462.
41. Elasfar A, Naeim H, Hussein M, et al. Percutaneous Retrieval of an Air Bullet From the Left Ventricle. *JACC Cardiovasc Interv.* 2018;11(13):1307-1308.
42. El-Damaty A, Love M, Parkash R. Detached tip of a transseptal sheath during left atrial ablation. *Catheter Cardiovasc Interv.* 2012;79(3):444-447.

43. Eleshra AS, Kolbel T, Rohlfes F, Scheerbaum M, Konstantinou N, Tsilimparis N. Emergent Use of a Branched Arch Device to Treat an Ascending Aortic Rupture: A Branch-to-Branch Through-and-Through Wire Technique to Compensate for Rotation Error. *J Endovasc Ther.* 2019;26(4):458-462.
44. Hennessey H, Valenti D, Cabrera T, Panet-Raymond V, Roberge D. Cardiac embolization of an implanted fiducial marker for hepatic stereotactic body radiotherapy: a case report. *Journal of Medical Case Reports.* 2009;3(1).
45. Jaber B, Arslan B, Ahmed O, Turba UC, Tasse J, Madassery S. Inadvertent noncaval filter deployment and its management. *J Vasc Surg Cases Innov Tech.* 2019;5(3):360-364.
46. Karim S, Abdelmessih M, Marieb M, Reiner E, Grubman E. Extraction of a Micra Transcatheter Pacing System: First-in-human experience. *HeartRhythm Case Rep.* 2016;2(1):60-62.
47. Morita J, Fukunaga M, Hiroshima K, Nagashima M, Miura M, Ando K. Retrieval of a Micra transcatheter pacing system in a heart with a preexisting lead. *Indian Pacing Electrophysiol J.* 2018;18(5):183-184.
48. Pasala TKR, Elmann EM, Ng AF, Song D, Ruiz CE. Transapical Retrieval of a Ruptured and Immobilized TAVR Balloon. *Semin Thorac Cardiovasc Surg.* 2019;31(3):470-472.
49. Sibona A, Smith JC, Srikureja DP, et al. Collaborative Management of Missile Injury to Right Ventricle and Subsequent Pulmonary Embolization. *Ann Thorac Surg.* 2018;106(6):e293-e294.
50. Vanezis AP, Wilson W, Smith W. Feasibility of 5 Fr contralateral radial access for hybrid chronic total occlusion percutaneous coronary intervention: Experience from a Single Tertiary Center in the United Kingdom. *Journal of Invasive Cardiology.* 2019;31(6):159-165.
51. Woodley-Cook J, Noel-Lamy M, Tan KT, Lindsay T. Use of an AMPLATZER duct occluder device to treat a recurrent arteriovenous fistula following endovenous laser treatment. *J Vasc Interv Radiol.* 2015;26(7):1082-1084.
52. Singla A, Foster M, Bennett WL, Young E, Attaran R, Mena-Hurtado C. Management of a patient with superior vena cava syndrome and a central venous catheter. *Vasc Dis Mgmt.* 2017;14(7):E157-E159.
53. Chiang WC, Ling-Lin Pai E, Chen PL, Shih CC, Chen IM. Endografting under Assistance of Transapical Body Floss Through-and-Through Wiring Technique and Rapid Ventricular Pacing for an Ascending Aortic Tuberculous Pseudoaneurysm. *Ann Vasc Surg.* 2017;45:270 e277-270 e211.
54. Hari P, Kirtane AJ, Bangalore S. Retrograde approach to an ostial left anterior descending chronic total occlusion through a left internal mammary artery graft. *Catheter Cardiovasc Interv.* 2016;87(6):E224-228.
55. Lawin D, Hansky B, Baumann P, Stellbrink C. Heart on a string: a novel approach to managing difficult access to the left pericardiacophrenic vein for phrenic nerve stimulation. *Herzschrittmacherther Elektrophysiol.* 2018;29(3):322-324.
56. Damluji AA, Kaynak HE, Heldman AW. Combined Retrograde/Antegrade Approach to Transcatheter Closure of an Aortic Paravalvular Leak. *Texas Heart Institute Journal.* 2015;42(5):443-447.
57. Sanghvi K, Selvaraj N, Luft U. Percutaneous closure of a perimembranous ventricular septal defect through arm approach (radial artery and basilic vein). *J Interv Cardiol.* 2014;27(2):199-203.
58. Beute TJ, Nolan MA, Merhi WM, Leung Wai Sang S. Use of EN Snare device for successful repositioning of the newest self-expanding transcatheter heart valve. *SAGE Open Med Case Rep.* 2018;6:2050313X18819933.

59. Kliger C, Al-Badri A, Wilson S, et al. Successful first-in-man percutaneous transapical-transseptal Melody mitral valve-in-ring implantation after complicated closure of a para-annular ring leak. *EuroIntervention*. 2014;10(8):968-974.
60. van Gelder BM, Houthuizen P, Bracke FA. Transseptal left ventricular endocardial pacing: preliminary experience from a femoral approach with subclavian pull-through. *Europace*. 2011;13(10):1454-1458.
61. Goodman SN, Berlin JA. The use of predicted confidence intervals when planning experiments and the misuse of power when interpreting results. *Ann Intern Med*. 1994;121:200-206.
62. AAMI/FDA Ad Hoc Risk Working G. *Postmarket risk management: A framework for incorporating benefit-risk assessments into correction and removal decisions*. Association for the Advancement of Medical Instrumentation;2016.
63. U.S. Food & Drug A. Factors to consider when making benefit-risk determinations in medical device premarket approval and de novo classifications: Guidance for industry and Food and Drug Administration staff. In: Administration USFD, ed. Rockville, MD: Center for Devices and Radiological Health; 2019.
64. Woodhouse JB, Uberoi R. Techniques for intravascular foreign body retrieval. *Cardiovasc Intervent Radiol*. 2013;36(4):888-897.
65. Kesselman AJ, Hoang NS, Sheu AY, Kuo WT. Endovascular Removal of Fractured Inferior Vena Cava Filter Fragments: 5-Year Registry Data with Prospective Outcomes on Retained Fragments. *J Vasc Interv Radiol*. 2018;29(6):758-764.
66. Johnstone RD, Stewart GA, Akoh JA, Fleet M, Akyol M, Moss JG. Percutaneous fibrin sleeve stripping of failing haemodialysis catheters. *Nephrol Dial Transplant*. 1999;14:688-691.
67. Reddy AS, Lang EV, Cutts J, Loh S, Rosen MP. Fibrin sheath removal from central venous catheters: an internal snare manoeuvre. *Nephrology Dialysis Transplantation*. 2007;22(6):1762-1765.
68. van Miert C, Hill R, Jones L. Interventions for restoring patency of occluded central venous catheter lumens. *Cochrane Database Syst Rev*. 2012(4):CD007119.
69. Paulus BM, Fischell TA. Retrieval devices and techniques for the extraction of intravascular foreign bodies in the coronary arteries. *J Interv Cardiol*. 2010;23(3):271-276.
70. Schechter MA, O'Brien PJ, Cox MW. Retrieval of iatrogenic intravascular foreign bodies. *J Vasc Surg*. 2013;57(1):276-281.
71. da Motta-Leal Filho JM, Carnevale FC, Nasser F, et al. Endovascular techniques and procedures, methods for removal of intravascular foreign bodies. *Rev Bras Cir Cardiovasc*. 2010;25(2):202-208.
72. Fazeny-Dorner B, Wenzel C, Berzlanovich A, et al. Central venous catheter pinch-off and fracture: recognition, prevention and management. *Bone Marrow Transplant*. 2003;31(10):927-930.
73. Hinke DH, Zandet-Stastny DA, Goodman LR, Quebbeman EJ, Krzywda EA, Andris DA. Pinch-off syndrome: A complication of implantable subclavian venous access devices. *Radiology*. 1990;177:353-356.
74. Altken DR, Minton JP. The "Pinch-Off Sign": A warning of impending problems with permanent subclavian catheters. *Am J Surg*. 1984;148:633-636.
75. Thompson KR. Nonsurgical retrieval of devices and foreign bodies. *Endovasc Today*. 2010;July:29-32.
76. Egglin TKP, Dickey KW, Rosenblatt M, Pollak JS. Retrieval of intravascular foreign bodies: Experience in 32 cases. *AJR*. 1995;164:1259-1264.

77. Wolf F, Schernthaner RE, Dirisamer A, et al. Endovascular management of lost or misplaced intravascular objects: experiences of 12 years. *Cardiovasc Intervent Radiol*. 2008;31(3):563-568.
78. Carroll MI, Ahanchi SS, Kim JH, Panneton JM. Endovascular foreign body retrieval. *J Vasc Surg*. 2013;57(2):459-463.
79. Torbicki A, Perrier A, Konstantinides S, et al. Guidelines on the diagnosis and management of acute pulmonary embolism: the Task Force for the Diagnosis and Management of Acute Pulmonary Embolism of the European Society of Cardiology (ESC). *Eur Heart J*. 2008;29(18):2276-2315.
80. Rentzos A, Lundqvist C, Karlsson JE, Vilmarsson V, Schnabel K, Wikholm G. Mechanical embolectomy for acute ischemic stroke in the anterior cerebral circulation: the Gothenburg experience during 2000-2011. *AJNR Am J Neuroradiol*. 2014;35(10):1936-1941.
81. Fanous AA, Siddiqui AH. Mechanical thrombectomy: Stent retrievers vs. aspiration catheters. *Cor et Vasa*. 2016;58(2):e193-e203.
82. Fisher RG, Ferreyro R. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies. *Am J Roentgenol*. 1978;130:541-548.
83. Cowley CG, Badran S, Rocchini AP, Lloyd TR. Snare-assisted vascular access: A new technique. *Cathet Cardiovasc Intervent*. 1999;47:315-318.
84. Kaufman JA, Kinney TB, Streiff MB, et al. Guidelines for the use of retrievable and convertible vena cava filters: report from the Society of Interventional Radiology multidisciplinary consensus conference. *J Vasc Interv Radiol*. 2006;17(3):449-459.
85. Sacks D, McClenny TE, Cardella JF, Lewis CA. Society of Interventional Radiology Clinical Practice Guidelines. *Journal of Vascular and Interventional Radiology*. 2003;14(9):S199-S202.
86. American College of R. American College of Radiology Appropriateness Criteria for the Management of Venous Thromboembolism and Inferior Vena Cava Filters. <https://acsearch.acr.org/docs/69342/Narrative/>. Published 2018. Accessed 6/27/2019, 2019.
87. National Kidney F. KDOQI clinical practice guideline for vascular access: 2018. https://www.kidney.org/sites/default/files/kdoqi_vasc-access-review2019_v2.pdf. Published 2018. Accessed 6/20/2019, 2019.

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