**COOK****MEDICAL****CE**  
0123**EN**  
**2****SteadySheath® Evolution® Tissue Stabilization Sheath**  
**SteadySheath® Evolution® Shortie Tissue Stabilization Sheath**  
Instructions for Use**CS**  
**6**Sheath pro stabilizaci tkáně **SteadySheath® Evolution®**  
Sheath pro stabilizaci tkáně **SteadySheath® Evolution® Shortie**  
Návod k použití**DA**  
**10****SteadySheath® Evolution® sheath til vævsstabilisering**  
**SteadySheath® Evolution® Shortie sheath til vævsstabilisering**  
Brugsanvisning**DE**  
**14****SteadySheath® Evolution® Gewebestabilisierungsschleuse**  
**SteadySheath® Evolution® Shortie**  
Gewebestabilisierungsschleuse  
Gebrauchsanweisung**EL**  
**19**Θηκάρι σταθεροποίησης ιστού **SteadySheath® Evolution®**  
Θηκάρι σταθεροποίησης ιστού **SteadySheath® Evolution®**  
**Shortie**  
Οδηγίες χρήσης**ES**  
**24**Vaina de estabilización de tejidos **SteadySheath® Evolution®**  
Vaina de estabilización de tejidos **SteadySheath® Evolution®**  
**Shortie**  
Instrucciones de uso**FR**  
**28**Gaine de stabilisation de tissus **SteadySheath® Evolution®**  
Gaine de stabilisation de tissus **SteadySheath® Evolution®**  
**Shortie**  
Mode d'emploi**HU**  
**32****SteadySheath® Evolution® szövetstabilizáló hüvely**  
**SteadySheath® Evolution® Shortie szövetstabilizáló hüvely**  
Használati utasítás**IT**  
**36**Guaina stabilizzatrice per tessuti **SteadySheath® Evolution®**  
Guaina stabilizzatrice per tessuti **SteadySheath® Evolution®**  
**Shortie**  
Istruzioni per l'uso**LT**  
**40**„**SteadySheath® Evolution®** audinių stabilizavimo vamzdelis  
„**SteadySheath® Evolution® Shortie**“ audinių  
stabilizavimo vamzdelis  
Naudojimo nurodymai**NL**  
**44****SteadySheath® Evolution® weefselstabiliserende sheath**  
**SteadySheath® Evolution® Shortie weefselstabiliserende sheath**  
Gebruiksaanwijzing**NO**  
**49****SteadySheath® Evolution® hylse til vevstabilisering**  
**SteadySheath® Evolution® Shortie hylse til vevstabilisering**  
Bruksanvisning**PL**  
**53**Koszulka **SteadySheath® Evolution®** do stabilizacji tkanki  
Koszulka **SteadySheath® Evolution® Shortie** do stabilizacji  
tkanki  
Instrukcja użycia**PT**  
**58**Bainha de estabilização tecidual **SteadySheath® Evolution®**  
Bainha de estabilização tecidual **SteadySheath® Evolution®**  
**Shortie**  
Instruções de utilização**RU**  
**62**Интердьюсер для стабилизации тканей  
**SteadySheath® Evolution®**  
Интердьюсер для стабилизации тканей  
**SteadySheath® Evolution® Shortie**  
Инструкция по применению**SV**  
**67****SteadySheath® Evolution® hylsa för vävnadsstabilisering**  
**SteadySheath® Evolution® Shortie hylsa för vävnadsstabilisering**  
Bruksanvisning

0 0 0 0 7 8 6 8 6 - 0 0 2

## SteadySheath® Evolution® Tissue Stabilization Sheath Instructions for Use

Federal law restricts this device to sale by or on the order of a physician.

### DEVICE DESCRIPTION

The SteadySheath Evolution Tissue Stabilization Sheath is comprised of a PTFE sheath with a stainless steel distal tip.

### INTENDED USE

The SteadySheath Evolution Tissue Stabilization Sheath is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects.

| Sheath Size (I.D.) | Compatible with Evolution or Evolution RL Sheath Sets | Compatible with Evolution Shortie or Evolution Shortie RL Sheath Sets |
|--------------------|-------------------------------------------------------|-----------------------------------------------------------------------|
| 9 Fr               | LR-TSS-9.0                                            | LR-TSS-SH-9.0                                                         |
| 11 Fr              | LR-TSS-11.0                                           | LR-TSS-SH-11.0                                                        |
| 13 Fr              | LR-TSS-13.0                                           |                                                                       |

### Use With Other Devices

The SteadySheath Evolution Tissue Stabilization Sheath may be used in conjunction with the following catheter/lead extraction devices from Cook:

- Locking Stylet
- Lead Extender
- Evolution Mechanical Dilator Sheath Set
- Evolution RL Controlled-Rotation Dilator Sheath Set
- Evolution Shortie Mechanical Dilator Sheath Set
- Evolution Shortie RL Controlled-Rotation Dilator Sheath Set

Be sure to closely follow the “Suggested Instructions for Use” for each device used.

### CONTRAINDICATIONS

None known

### WARNINGS

- **When using dilator sheaths or sheath sets do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.**
- **When using a Locking Stylet:**
  - **Do not abandon a catheter/lead in a patient with a Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.**
  - **Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.**
  - **Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.**
- **Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:**
  - **The item to be removed is of a dangerous shape or configuration,**
  - **The likelihood of catheter/lead disintegration resulting in fragment embolism is high or vegetations are attached directly to the catheter/lead body.**
- **Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.**
- **Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.**

- **When advancing sheaths including the Evolution device and/or SteadySheath Tissue Stabilization Sheath, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to avoid damage to vessel walls.**
- **Excessive force with sheaths (including the Evolution device and/or SteadySheath Tissue Stabilization Sheath) used intravascularly may result in damage to the vascular system requiring surgical repair.**
- **Evolution Shortie or Evolution Shortie RL should only be used to minimally enter the vessel. Do not attempt to negotiate the Evolution Shortie or Evolution Shortie RL device past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur.**

## **PRECAUTIONS**

- Prior to the procedure, consider the size of the catheter/lead in relation to the size of the Lead Extraction™ devices to determine possible incompatibility.
- If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the non-targeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following:

## **PREPARATIONS**

- Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.
- Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.
- Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position. Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray.
- Echocardiography and cardiothoracic surgery should be rapidly available.
- Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.
- Establish back up pacing as necessary.
- Have available an extensive collection of sheaths, Lead Control Devices (Locking Stylet and Lead Extender), stylets to unscrew active fixation leads, snares, and accessory equipment.

## **PROCEDURE**

- Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.
- If using sheaths or sheath sets, including the Evolution device and/or the SteadySheath Tissue Stabilization Sheath, the following precautions should be followed:
  - Prior to using sheaths including the Evolution device or SteadySheath it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.
  - The Evolution Shortie, Evolution Shortie RL or Byrd Telescoping Stainless Steel Dilator Sheath Set should be used only to minimally enter the vessel.
  - Be aware that insertion of a Byrd Telescoping Stainless Steel Dilator Sheath Set over a Byrd Dilator Sheath Set Polypropylene or the Evolution device may damage the inner sheath.
  - If excessive scar tissue or calcification prevents safe advancement of sheaths, consider an alternate approach.
  - If the catheter/lead breaks, evaluate fragment; retrieve as indicated.
  - If hypotension develops, rapidly evaluate; treat as appropriate.
- Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

## POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

dislodging or damaging nontargeted catheter/lead  
chest wall hematoma  
thrombosis  
arrhythmias  
acute bacteremia  
acute hypotension  
pneumothorax  
stroke  
migrating fragment from catheter/object  
pulmonary embolism  
laceration or tearing of vascular structures or the myocardium  
hemopericardium/pericardial effusion  
cardiac tamponade  
hemothorax  
cardiac arrest  
death

## INSTRUCTIONS FOR USE

### Suggested Instructions for Use

#### SteadySheath Evolution Tissue Stabilization Sheath

1. Surgically expose the proximal end of the indwelling catheter/lead and remove the catheter/lead from its connections (if connected). Remove all suture and tie-down materials.
2. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.
3. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature or affix a compression coil (Cook) at the proximal end of the catheter/lead, compressing the insulation against the coil and Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle or to the suture tie loop.  
  
NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.
4. For an active fixation catheter/lead, attempt to unscrew the catheter/lead by rotating the catheter/lead and Locking Stylet counterclockwise if appropriate.
5. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Locking Stylet and catheter/lead to remove it.
6. If the catheter/lead is not removed from the vessel with gentle pulling then using dilator sheaths, including the Evolution device, Evolution Shortie device, SteadySheath or other retrieval devices may help separate the catheter/lead from tissue encapsulation as follows: Advance dilator sheaths (or other retrieval device) over the proximal wire extender of the Locking Stylet (if applicable).

7. With the SteadySheath placed outside the appropriate inner sheath for telescopic action, insert the proximal free end of the catheter/lead with Locking Stylet in place, into the distal end of the inner Evolution device or Evolution Shortie device. Advance the catheter/lead until it completely exits the opposite (proximal) end of the sheath set.

8. Apply adequate retracting pressure or tension on the catheter/lead and/or its wire guide or Locking Stylet (if used). This is critical to safe passage of the sheath set over the catheter/lead. If tension is inadequate, the catheter/lead may buckle, precluding the sheath set from advancement along the appropriate path.

9. With the catheter/lead in tension, and under fluoroscopic guidance, advance the inner sheath along the length of the catheter/lead and, if applicable, into the blood vessel. Advance the SteadySheath over the inner sheath keeping it within the vessel.

**WARNING: Evolution Shortie or Evolution Shortie RL should only be used to minimally enter the vessel. Do not attempt to negotiate the Evolution Shortie or Evolution Shortie RL device past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur.**

10. To utilize SteadySheath, hold it stationary while advancing and/or rotating the inner sheath to avoid undue tissue movement. This technique will aid in preventing adjacent tissue and indwelling catheters from rotating around the inner sheath when it is advancing.

11. For cardiac leads, if the lead has not been freed by the time the sheaths near the myocardium, position the SteadySheath so that the blunt end is directed toward the myocardium. Pull the inner sheath back several centimeters.

Apply countertraction with the SteadySheath as follows:


Firmly hold the outer sheath one centimeter from the heart wall and steadily pull the Locking Stylet back; the pacing tip will be pulled into the sheath.

12. When the catheter/lead has become unentangled and freed from tissue attachment, carefully remove from the patient, the catheter/lead, Locking Stylet if used and SteadySheath (or other retrieval device).

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Vandergrift, PA 15690 U.S.A.
COOK MEDICAL EUROPE LTD