- EN 2 **Needle's Eye Snare®**
Instructions for Use
- CS 7 **Očko Needle's Eye Snare®**
Návod k použití
- DA 12 **Needle's Eye Snare® (slyng)**
Brugsanvisning
- DE 17 **Needle's Eye Snare® Schlinge**
Gebrauchsanweisung
- EL 23 **Βρόχος Needle's Eye Snare®**
Οδηγίες χρήσης
- ES 29 **Asa Needle's Eye Snare®**
Instrucciones de uso
- FR 35 **Récupérateur Needle's Eye Snare®**
Mode d'emploi
- HU 41 **Needle's Eye Snare®**
Használati utasítás
- IT 47 **Needle's Eye Snare®**
Istruzioni per l'uso
- NL 53 **Needle's Eye Snare®**
Gebruiksaanwijzing
- PL 59 **Pętla Needle's Eye Snare®**
Instrukcja użycia
- PT 65 **Laço Needle's Eye Snare®**
Instruções de utilização
- RU 71 **Петля Needle's Eye Snare®**
Инструкция по применению
- SV 77 **Needle's Eye Snare®-snaran**
Bruksanvisning

Needle's Eye Snare®

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Needle's Eye Snare® is comprised of a stainless steel/nitinol wire retrieval mechanism housed within two telescoping TFE sheaths, a protective 8 Fr sheath and a removable 12 Fr straight inner femoral sheath that may be advanced forward to close the retrieval mechanism. The retrieval mechanism consists of a "needle's eye" loop and "threader." The snare is used through a 16 Fr TFE introducer sheath which incorporates a hemostasis valve into its proximal end and has either a straight or curved distal segment. Additional introducer set components (needle, wire guide, dilator) are used for vessel entry.

For model LR-NES541 only, a length appropriate 12 Fr Curved Inner Sheath is included.

For models: LR-NES001, LR-NES002, a 12 Fr Curved Inner Femoral Sheath and a Curved 16 Fr Introducer Sheath Set are available separately.

Catalog Number	Product
LR-NES001, LR-NES002	Needle's Eye Snare Retrieval Sets
LR-NES541	Needle's Eye Snare Retrieval Set - 54 cm
LR-CRVFEM001	12 Fr Curved Inner Femoral Sheath
LR-CSS16	Curved 16 Fr Introducer Sheath Set - Femoral

INDICATIONS

The Needle's Eye Snare Retrieval Sets, 12 Fr Curved Inner Femoral Sheath, Curved 16 Fr Introducer Sheath Set Femoral are intended for use in patients requiring the percutaneous retrieval of indwelling catheters, cardiac leads, fragments of catheter tubing or wire guides, and other foreign objects.

CONTRAINDICATIONS

None known

WARNINGS

Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:

- the item to be removed is of a dangerous shape or configuration,
- the likelihood of catheter/lead disintegration resulting in fragment embolism is high,
- or vegetations are attached directly to the catheter/lead body.

Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.

Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.

PRECAUTIONS

Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility. If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the nontargeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

Do not use excessive force to grasp and pull large portions of doubled-over lead/catheter into the 12 Fr or 16 Fr sheath, as it could become lodged, preventing further movement of the sheaths or snare(s).

Avoid using the Needle's Eye Snare in a manner that could potentially put unwarranted forces and fatigue on the wire. The device is not intended to be excessively twisted and bent while using. Such techniques may cause integrity issues leading to wire breakage within the patient.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following precautions:

Preparations

- Obtain a thorough patient history, including patient blood type.
- Appropriate blood products should be rapidly available.

Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed. Perform radiographic/ echocardiographic evaluation of catheter/lead condition, type, and position.

Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray. Echocardiography and cardiothoracic surgery should be rapidly available. Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach.

Establish back-up pacing as necessary.

Have available an extensive collection of sheaths, Locking Stylets, stylets to unscrew active fixation leads, snares, and accessory equipment.

Procedure

Use fluoroscopic monitoring during ALL sheath, catheter/lead, or snare manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.

Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

When advancing sheath(s) over the catheter/lead, **maintain adequate tension on the catheter/lead** (via the basket snare or Needle's Eye Snare®) to guide the sheath(s) and prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of sheath(s), consider an alternate approach.

Excessive force with sheath(s) used intravascularly may result in damage to the vascular system requiring surgical repair.

If the catheter/lead breaks, evaluate fragment; retrieve as indicated.

If hypotension develops, rapidly evaluate; treat as appropriate.

Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

ADVERSE EVENTS

Potential Adverse Events

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

- dislodging or damaging nontargeted catheter/lead
- chest wall hematoma
- thrombosis
- arrhythmias
- acute bacteremia
- acute hypotension
- pneumothorax
- stroke
- migrating fragment from catheter/object
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium/pericardial effusion
- cardiac tamponade
- hemothorax
- cardiac arrest
- death

INSTRUCTIONS FOR USE

Suggested Instructions for Use: Needle's Eye Snare®

NOTE: In addition to closely following the "Suggested Instructions for Use" for this device, refer to the "Clinical Considerations for Removing Leads, Catheters, or Other Indwelling Objects" shown in the next section of this insert.



NOTE: Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

CAUTION: Use fluoroscopic monitoring during all sheath, catheter/lead, or snare manipulations.

CAUTION: Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

1. Use the Seldinger technique for sterile insertion of a wire guide into the vessel. Under fluoroscopic guidance, insert the 16 Fr introducer sheath.
2. Review and complete this step only if replacement of the existing 12 Fr inner sheath with a new straight or curved 12 Fr inner sheath is required. Otherwise, go to step 3.

a) On the proximal base of the Needle's Eye Snare® assembly, release the introducer tube from its storage position by rotating the Luer connector.


b) Slide the introducer tube forward, toward the distal tip of the Needle's Eye Snare®.



Introducer Tube Advanced

c) Carefully advance the introducer tube over the Needle's Eye Snare® loop, thereby drawing the loop inside the introducer tube. Leave a short section of the loop exposed.

Loop Drawn Inside Introducer Tube

d) Insert the distal end of the introducer tube into the proximal Luer connector of the 12 Fr Inner Femoral Sheath.


Introducer Tube Inserted into Luer Connector

e) Push the Needle's Eye Snare® assembly through the introducer tube and into the 12 Fr inner sheath. Advance the Needle's Eye Snare® assembly until the tip of the loop begins to emerge from the 12 Fr inner sheath.

f) Return the introducer tube to its storage position at the proximal base of the Needle's Eye Snare® assembly, and secure its position by rotating the Luer connector.

3. Introduce the Needle's Eye Snare® assembly (which is now retracted within the 12 Fr inner sheath) into the 16 Fr introducer sheath already in place in the vessel. Advance the assembly, using fluoroscopic guidance, until the distal end of the 12 Fr inner sheath, with Needle's Eye Snare®, has exited the distal end of the 16 Fr introducer sheath and is adjacent to the indwelling object. The side arm flush port may be used for fluid infusion or to prefill the snare's inner lumen with saline.

As the sheath is retracted the "needle's eye" loop begins to emerge, and curl into its preformed shape.



Once positioned at site of foreign object, the "needle's eye" is expelled from the sheath.

4. Deploy the "needle's eye," while using fluoroscopic monitoring, by retracting the 12 Fr inner sheath. As the sheath is retracted, the "needle's eye" will emerge and begin to curl into its preformed shape. Ensure that the "needle's eye" loop positions itself around the indwelling object.


"Needle's Eye"
"Threader"


5. Once the "needle's eye" loop is positioned over the body of the indwelling object, the "threader" should be fully extended by pushing the black knob on the handle completely forward. Ensure that the body of the indwelling device is totally encircled by the "needle's eye" loop and "threader." **Engagement occurs when the "threader" is fully positioned inside the loop of the "needle's eye," with the black slide knob of the handle fully forward.**

6. Once engagement between the "needle's eye" and "threader" is confirmed, the indwelling object can be firmly grasped or captured by holding the handle of the Needle's Eye Snare® stationary and advancing the 12 Fr sheath distally until the indwelling object is held firmly against the distal end of the 12 Fr sheath.

7. The indwelling object is prepared for removal by advancing the 16 Fr introducer sheath over the object and/or retracting the 12 Fr inner sheath, until the object enters the 16 Fr sheath.

8. When the indwelling object has been freed and retracted into the 16 Fr introducer sheath, withdraw the Needle's Eye Snare® assembly with attached object out of the vessel, while maintaining hemostasis.

Clinical Considerations for Removing Leads, Catheters or Other Indwelling Objects

Clinical experience in the removal of leads has identified several considerations for lead removal techniques using superior or femoral approaches. Physicians highly experienced with lead removal techniques have suggested the following considerations for removal of leads, catheters, or other indwelling objects.

Part I. Accessing the Vessel

CAUTION: Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

CAUTION: Use fluoroscopic monitoring during all sheath, catheter/lead, or snare manipulations.

CAUTION: Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

1. Insert the dilator into the 16 Fr sheath (outer sheath). Use the Seldinger technique for sterile insertion of the syringe-mounted needle and wire guide into the vein. Being certain that the dilator leads the sheath, insert the 16 Fr introducer sheath and dilator under fluoroscopic guidance. Remove the dilator and advance the sheath until it extends to the middle third of the vena cava.

Part II. Snaring the Indwelling Object

2. Follow the "Suggested Instructions for Use" to:
advance the 12 Fr sheath containing the snare to be used (Tip Deflecting Wire Guide with Dotter Basket, or the Needle's Eye Snare®) inside the 16 Fr introducer sheath, deploy the snare, and snare the catheter/lead or other indwelling object.

Part III. Removing the Indwelling Object

3. If removing an active fixation lead, rotate the snare (helical loop basket or Needle's Eye Snare®) counterclockwise to unscrew the fixation mechanism.
4. While maintaining adequate tension on the indwelling object via the closed snare (basket or Needle's Eye Snare®), prepare the indwelling object for removal by advancing the 16 Fr sheath over the object to encapsulate it. Cautiously disrupt fibrotic adhesions encountered during sheath advancement. (Where possible, the 12 Fr sheath can also be advanced to aid in dilation.)
5. For a cardiac lead, if the lead has not been freed by the time the sheath nears the myocardium, position the 16 Fr sheath about one centimeter from the myocardium, making sure that the 12 Fr sheath has been withdrawn into the 16 Fr sheath at least 6 cm. As the 16 Fr sheath is held in position by stabilizing the sheath at its venous entry site in the groin, apply countertraction by adding tension to the snare (basket or Needle's Eye Snare®), thereby pulling the lead tip into the sheath and drawing the endocardium against the 16 Fr sheath's distal tip. Steadily pull on the lead until the lead tip is freed from the myocardium and is pulled into the sheath.
6. When the indwelling object has been freed from tissue entrapment and retracted into the 16 Fr sheath, withdraw the sheath back into the vena cava. Remove the entire ensemble while maintaining hemostasis.

NOTE: If removing multiple catheters/leads, leave the 16 Fr sheath in place until all catheters/leads have been extracted. If using the Tip Deflecting Wire Guide and Dotter Basket, a new Tip Deflecting Wire Guide, Dotter Basket, and 12 Fr sheath should be used to extract subsequent leads.

NOTE: To prevent potential complications, carefully examine any broken catheters/leads. A long exposed lead coil could lacerate the heart, or a loose fragment could migrate, possibly into the pulmonary artery. A pigtail catheter or closed loop snare may be used to retrieve a catheter/lead from the pulmonary artery, repositioning it in the ventricle for removal.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2021
Vandergrift, PA 15690 U.S.A.
COOK MEDICAL EUROPE LTD