



Protection through Innovation

From waste management to sharps protection


to medication error prevention, customers rely

on Merit Medical when safety matters most.



MERIT SAFETY SOLUTIONS

Merit Medical is committed to providing unique medical devices that will assist clinical personnel in their efforts to increase safety.



Medication Error Prevention



PAL™ (PEN AND LABELS)

The PAL is a streamlined medication labeling system. Each label is made of transparent moisture resistant material. The sterile smear resistant permanent marker is provided in order to print pertinent information.



MEDALLION® SYRINGES

These color-coded polycarbonate syringes help to minimize the potential risk of medication errors. The clear barrel provides easy identification of air bubbles with a "smooth as glass" feel.

Sharps Protection



A•S•K MERIT ACCESS SAFETY KIT®

The Merit access safety kit incorporates shielded technology that helps minimize the risk of sharps injuries associated with gaining vascular access. The A•S•K features the SecureLoc introducer needle and the Futura safety scalpel.



FUTURA® SAFETY SCALPELS

These automated retractable safety scalpels are designed to protect clinicians and patients from accidental cuts. The one-handed activation allows the clinician to focus on the surgical site.



SECURELOC™ ANGIOGRAPHIC INTRODUCER NEEDLE

The unique design of the SecureLoc introducer needle allows clinicians to easily activate a low profile safety mechanism, shielding the needle after use to minimize any harmful needlesticks.



SHORTSTOP® TEMPORARY SHARPS HOLDER

The ShortStop temporarily holds sharps in a protective environment between multiple sticks during a procedure. Designed with colorful bullseye images and an adhesive strip for the securement to a drape, cover or tray.

SHORTSTOP ADVANTAGE™ TEMPORARY SHARPS HOLDER


The ShortStop Advantage has a new side access feature allowing placement of a loaded syringe attached to a needle thus preventing any tipping or needle bending.

Waste Management



BACKSTOP® & BACKSTOP® PLUS CLOSED WASTE BASINS

The OSHA compliant BackStop and BackStop Plus are designed with particulate-free materials, a seal tight lid and a 500ml fluid capacity. Designed for quick and safe drainage elimination.



MINISTOP® & MINISTOP® PLUS CLOSED WASTE BASINS

The OSHA compliant MiniStop and MiniStop Plus are designed with particulate-free materials, a seal tight lid and a 100ml fluid capacity. Designed for smaller drainage procedures.



DUGOUT® CLOSED WASTE BASIN

The OSHA compliant DugOut is designed with particulate-free material, a secondary compartment for storage of ancillary products or saline and a seal tight lid with a 1000ml fluid capacity. Designed for quick and safe elimination of the largest volume drainage procedures.



TRIPLEPLAY™ CLOSED WASTE BASIN

This OSHA compliant disposal system is designed to reduce the exposure to harmful bloodborne pathogens. The TriplePlay comes equipped with a covered waste basin, an absorbent catch basin and an ancillary product tray.



MERIT® DISPOSAL DEPOT

The Merit closed waste management system is a cost-effective, simple way to satisfy OSHA guidelines for managing blood waste.



CAPTIVA® BLOOD CONTAINMENT DEVICE

Attached to an angiographic needle, the Captiva minimizes bloodborne hazards associated with arterial needle sticks. The design allows clinicians to visualize blood flashback while limiting their exposure to harmful bloodborne pathogens.

**SHORTSTOP ADVANTAGE
TEMPORARY SHARPS HOLDER DEVICE**

SS20 - Temporary Sharps Holder Device
25 Units Per Box

SHORTSTOP TEMPORARY SHARPS HOLDER DEVICE

SS10 - Temporary Sharps Holder Device
25 Units Per Box

FUTURA SAFETY SCALPELS

SMS210 - #10 Blade
SMS211 - #11 Blade
SMS215 - #15 Blade
50 Units Per Box

SECURELOC ANGIOGRAPHIC INTRODUCER NEEDLE

SL18T71W - 18G, 7cm (2 3/4")
SL18T71WS - 18G, 7cm (2 3/4") Seldinger Shield
SL19T71W - 19G, 7cm (2 3/4")
SL19T71WS - 19G, 7cm (2 3/4") Seldinger Shield
SL21T71W - 21G, 7cm (2 3/4")
SL21T71WS - 21G, 7cm (2 3/4") Seldinger Shield
25 Units Per Box

A•S•K MERIT ACCESS SAFETY KIT

ASK200 - (featuring the #11 Futura & 18G, 7cm SecureLoc)
25 Units Per Box

MEDALLION SYRINGES

Seven colors available in 1ml, 3ml, 6ml, 10ml, 20ml, 30ml and 60ml product configurations. Custom printed syringes available by request.

See product catalog for additional information.
25 Units Per Box

PAL (PEN AND LABEL)

PAL100 - Pen and 35 Label Set
PAL200 - Pen and 24 Label Set
PAL225 - Pen, 24 Label Set & 5 Guide Wire Bowl Label Set
PAL250 - Pen, 24 Label Set & 7 Large General Bowl Label Set
PAL300 - Pen and 12 Label Set
PAL325 - Pen, 12 Label Set & 5 Guide Wire Bowl Set
PAL350 - Pen, 12 Label Set & 7 Large General Bowl Label Set
PAL400 - Pen and 12 Label Set (Blank)
PAL425 - Pen, 12 (Blank) Label Set & 5 Guide Wire Bowl Set
PAL450 - Pen, 12 (Blank) Label Set & 7 Large General Bowl Set
PAL500 - 5 Guide Wire Bowl Set
PAL600 - 7 Large General Bowl Set
PAL700 - Pen, 35 Label Set & 7 Large General Bowl Set - NO INITIALS & DATE

Custom PAL label sets available by request.
25 Units Per Box



MERIT DISPOSAL DEPOTS

MDD100 - PVC bag, 72" (183cm) large bore tubing with fixed male luer.

MDD110 - PVC bag, 72" (183cm) large bore tubing, drape clip and one-way check valve.

MDD200 - PVC bag, 48" (122cm) large bore tubing. Includes check relief valve with female luer and fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented spike.

MDD220 - PVC bag, 48" (122cm) large bore tubing. Includes check relief valve with female luer and fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented spike, with optional reflux valve.

MDD300 - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve with rotating luer. Fluid administration set with 72" (183cm) large bore tubing, roller clamp and user-selectable vented spike.

MDD400 - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve with rotating luer. Fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented drip chamber with micro filter.

MDD500 - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve, rotating luer and standard pressure 3-way Marquis Stopcock. Two fluid administration sets with 72" (183cm) large bore tubing, roller clamp and user selectable vented spikes.

MDD520 - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve, rotating luer and standard pressure 3-way Marquis Stopcock. Two fluid administration sets with 72" (183cm) large bore tubing, roller clamp and user selectable vented spikes, with optional reflux valve.

Custom configurations available upon request.
20 Units Per Box

DUGOUT 1000ML CLOSED WASTE BASIN

DUG100 - Closed Waste Basin
DUG1SS - Closed Waste Basin (contains the SS10)
6 Units Per Box

BACKSTOP & BACKSTOP PLUS 500ML CLOSED WASTE BASINS

MBS100 - Closed Waste Basin
MBSP100 - Closed Waste Basin with Temporary Sharps Holder
5 Units Per Box

MINISTOP & MINISTOP PLUS 100ML CLOSED WASTE BASINS

MMS100 - Closed Waste Basin
MMSP100 - Closed Waste Basin with Temporary Sharps Holder
20 Units Per Box

TRIPLEPLAY 200ML CLOSED WASTE BASIN

MTP100 - Closed Waste Basin
10 Units Per Box

CAPTIVA BLOOD CONTAINMENT DEVICE

CAN100 - Blood Containment Device
25 Units Per Box

For additional information, consult your Merit Medical sales representative or visit our website at www.merit.com.

All products are latex free.