



Catch the Vision

MERT **MAKIN**




MINI ACCESS KIT

*Merit strikes “gold” again with the MAKNV.
Excellent transitions for a smooth introduction.*


ORDERING INFORMATION

Catalog Number Contents

- MAK-NV-001** 0.038" (0.97mm) X 150cm double ended heavy duty 3mm J guide wire
0.018" (0.46mm) X 60cm stainless steel guide wire with platinum tip
6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
21g X 15cm introducer needle with trocar stylet
- MAK-NV-002** 0.038" (0.97mm) X 150cm double ended heavy duty 3mm J guide wire
0.018" (0.46mm) X 60cm nitinol guide wire with platinum tip
6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
21g X 15cm introducer needle with trocar stylet
- MAK-NV-003** 0.018" (0.46mm) X 60cm stainless steel guide wire with platinum tip
6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
21g X 15cm introducer needle with trocar stylet
- MAK-NV-004** 0.018" (0.46mm) X 60cm nitinol guide wire with platinum tip
6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
21g X 15cm introducer needle with trocar stylet
- MAK-NV-005*** 21g X 15cm introducer needle with trocar stylet
- MAK-NV-006*** 6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
21g X 15cm introducer needle with trocar stylet
- MAK-NV-007*** 6F (2.0mm) X 20cm co-axial introducer with dilator and stiffening cannula
- VA-40** Valve Adapter - accepts up to a 4F device
- MAK005-60** 0.018" (0.46mm) X 60cm stainless steel guide wire with platinum tip
- MAK001N-60** 0.018" (0.46mm) X 60cm nitinol guide wire with platinum tip
Packaged 10 per box * Packaged 5 per box



VALVE ADAPTER



REVOLUTION®



RESOLVE® LOCKING



DRAINAGE DEPOT™



RESOLVE® BILIARY

FEATURES & BENEFITS

- RADIOPAQUE MARKER TIP**
clearly shows tip of introducer
- SUPERIOR TIP DESIGN**
decreases buckling and introducer peel back during insertion
- ECHO ENHANCED NEEDLE**
increases visualization for ultrasound placement
- 4F INTRODUCER ID**
allows 0.018" guide wire may be left in place for safety during larger wire placement

MERIT'S DRAINAGE FAMILY

offers all you need from introducers to drainage catheters and fixation devices to bags and accessories.

