




Protection through Innovation

*From waste management to sharps protection
to medication error prevention, customers rely
on Merit Medical when safety matters most.*

MERIT SAFETY SOLUTIONS

Merit Medical is committed to providing unique medical devices that will assist clinical personnel in their efforts to increase safety.



Medication Error Prevention



PAL™ (PEN AND LABELS)

The PAL is a streamlined medication labeling system. Each label is made of transparent moisture resistant material. The sterile smear resistant permanent marker is provided in order to print pertinent information.



MEDALLION® SYRINGES

These polycarbonate syringes help to minimize the potential risk of medication errors. The clear barrel provides easy identification of air bubbles with a "smooth as glass" feel. Select from options of frosted field, colored plungers, preprinted and custom clichés.

Sharps Protection



FUTURA® SAFETY SCALPELS

These automated retractable safety scalpels are designed to protect clinicians and patients from accidental cuts. The one-handed activation allows the clinician to focus on the surgical site.

Waste Management



BACKSTOP® & BACKSTOP® PLUS CLOSED WASTE BASINS

The OSHA compliant BackStop and Back-Stop Plus are designed with particulate-free materials, a seal tight lid and a 500ml fluid capacity. Designed for quick and safe fluid elimination.



SECURELOC™ ANGIOGRAPHIC

Introducer Needle

The unique design of the SecureLoc introducer needle allows clinicians to easily activate a low profile safety mechanism, shielding the needle after use to minimize any harmful needlesticks.



MINISTOP® & MINISTOP® PLUS CLOSED WASTE BASINS

The OSHA compliant MiniStop and MiniStop Plus are designed with particulate-free materials, a seal tight lid and a 100ml fluid capacity. Designed for smaller fluid procedures.



SHORTSTOP® TEMPORARY SHARPS HOLDER

The ShortStop temporarily holds sharps in a protective environment between multiple sticks during a procedure. Designed with colorful bullseye images and an adhesive strip for the securement to a drape, cover or tray.



DUGOUT® CLOSED WASTE BASIN

The OSHA compliant DugOut is designed with particulate-free material, a secondary compartment for storage of ancillary products or saline and a seal tight lid with a 1000ml fluid capacity. Designed for quick and safe elimination of the largest volume drainage procedures.



SHORTSTOP ADVANTAGE®

The ShortStop Advantage has a new side access feature allowing placement of a loaded syringe attached to a needle thus preventing any tipping or needle bending.



MERIT® DISPOSAL DEPOT

The Merit closed waste management system is a cost-effective, simple way to satisfy OSHA guidelines for managing blood waste.

SHORTSTOP GRANDSTAND™

The best seat in the house that holds your sharp devices upright while utilizing minimal space.

SHORTSTOP TEMPORARY SHARPS HOLDER DEVICES

SS10 - ShortStop -25 Units Per Box

SS20 - ShortStop Advantage - 25 Units Per Box

SS30 - ShortStop Grandstand - 10 Units Per Box

FUTURA SAFETY SCALPELS

SMS210 - #10 Blade

SMS211 - #11 Blade

SMS215 - #15 Blade

50 Units Per Box

SECURELOC ANGIOGRAPHIC INTRODUCER NEEDLE

SL18T71W - 18G, 7cm (2 3/4")

SL18T71WS - 18G, 7cm (2 3/4") Seldinger Shield

SL21T71W - 21G, 7cm (2 3/4")

SL21T71WS - 21G, 7cm (2 3/4") Seldinger Shield

25 Units Per Box

MEDALLION SYRINGES

MSSF011 - 1ml Frosted Field Syringe

MSSF031 - 3ml Frosted Field Syringe

MSSF061 - 6ml Frosted Field Syringe

MSSF111 - 10ml Frosted Field Syringe

MSSF121 - 20ml Frosted Field Syringe

MSSF131 - 30ml Frosted Field Syringe

MSSF160 - 60ml Frosted Field Syringe

Seven colors available in 1ml, 3ml, 6ml, 10ml, 20ml, 30ml and 60ml product configurations. Custom printed syringes available by request.

See product catalog for additional information. - 25 Units Per Box

PAL (PEN AND LABEL)

PAL100 - Pen and 35 Label Set

PAL200 - Pen and 24 Label Set

PAL225 - Pen, 24 Label Set & 5 Guide Wire Bowl Label Set

PAL250 - Pen, 24 Label Set & 7 Large General Bowl Label Set

PAL300 - Pen and 12 Label Set

PAL325 - Pen, 12 Label Set & 5 Guide Wire Bowl Set

PAL350 - Pen, 12 Label Set & 7 Large General Bowl Label Set

PAL400 - Pen and 12 Label Set (Blank)

PAL425 - Pen, 12 (Blank) Label Set & 5 Guide Wire Bowl Set

PAL450 - Pen, 12 (Blank) Label Set & 7 Large General Bowl Set

PAL500 - 5 Guide Wire Bowl Set

PAL600 - 7 Large General Bowl Set

PAL700 - Pen, 35 Label Set & 7 Large General Bowl Set - NO INITIALS & DATE

Custom PAL label sets available by request. - 25 Units Per Box

MERIT DISPOSAL DEPOTS

MDD100 - PVC bag, 72" (183cm) large bore tubing with fixed male luer.

MDD110 - PVC bag, 72" (183cm) large bore tubing, drape clip and one-way check valve.

MDD200H - PVC bag, 48" (122cm) large bore tubing. Includes check relief valve with female luer and fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented spike.

MDD220H - PVC bag, 48" (122cm) large bore tubing. Includes check relief valve with female luer and fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented spike, with optional reflux valve.

MDD300H - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve with rotating luer. Fluid administration set with 72" (183cm) large bore tubing, roller clamp and user-selectable vented spike.

MDD400H - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve with rotating luer. Fluid administration set with 72" (183cm) large bore tubing, roller clamp and vented drip chamber with micro filter.

MDD500H - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve, rotating luer and standard pressure 3-way Marquis Stopcock. Two fluid administration sets with 72" (183cm) large bore tubing, roller clamp and user selectable vented spikes.

MDD520H - PVC bag and 48" (122cm) large bore tubing. Includes check relief valve, rotating luer and standard pressure 3-way Marquis Stopcock. Two fluid administration sets with 72" (183cm) large bore tubing, roller clamp and user selectable vented spikes, with optional reflux valve.

Custom configurations available upon request.

20 Units Per Box

DUGOUT 1000ML CLOSED WASTE BASIN

DUG100 - Closed Waste Basin

DUG1SS - Closed Waste Basin (contains the SS10)

6 Units Per Box

BACKSTOP & BACKSTOP PLUS 500ML CLOSED WASTE BASINS

MBS100 - Closed Waste Basin

MBSP100 - Closed Waste Basin with Temporary Sharps Holder

5 Units Per Box

MINISTOP & MINISTOP PLUS 100ML CLOSED WASTE BASINS

MMS100 - Closed Waste Basin

MMSP100 - Closed Waste Basin with Temporary Sharps Holder

20 Units Per Box

For additional information, consult your Merit Medical sales representative or visit our website at www.merit.com.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

400848001/B ID 052410