

MERTMEDICAL[®]

| dialysis


contents


HEMODIALYSIS

prepare
access
navigate
implant
treat
integrate

PERITONEAL DIALYSIS

prepare
place
select
integrate
assist
complete





Merit Medical's full suite of innovative dialysis products can help you improve your patient's quality of life when treating end-stage renal disease (ESRD). This comprehensive line of peritoneal dialysis and hemodialysis catheters, implantation kits, and accessories has been designed to meet the unique clinical needs of dialysis treatment.

hemodialysis

prepare


access


navigate

implant

treat

integrate







ShortStop Advantage®
Temporary Sharps Holder



Medallion®
Syringe

PAL™
Medication
Labeling System



➤ promote safe preparation

Merit Medical has a complete line of products designed to promote the safety of clinicians and patients. **SHIELDED** scalpels, needles, and **SECURE** temporary sharps holders help protect clinicians from accidental sharps injuries. Multiple color options and smear-resistant labeling for syringes and containers help you prevent medication errors.



Futura®
Safety Scalpel

FlowGuard™ *
Valved Introducer

Merit Advance®
Straight & Curved Needle

access

➤ start with quality

Achieve **QUICK** and effective access with our **SHARP** angiographic needles, scalpels, mini access kits designed to lessen trauma and reduce bleeding, and sheath introducers with excellent transitions for **SMOOTH** insertion.


**FlowGuard is a trademark of Greatbatch, Inc.*

navigate

ProGuide™
Tri-Ball Tunneler

Impress®
30 cm Diagnostic Peripheral Catheter

CentrosFLO®
Tunneler





InQwire®
Diagnostic Guide Wire




Merit Laureate®
Hydrophilic Guide Wire



navigate with confidence

LUBRICIOUS, hydrophilic guide wires reduce friction and enhance your ability to make rapid catheter exchanges. The **TAPERED** tips, kink resistant hub design, and unique braided shafts give Merit's diagnostic catheters both flexibility and **STEERABILITY**. Our diagnostic guide wires are manufactured with a unique PTFE precoating process which results in a consistently smooth surface for smoother tracking through difficult vasculature.

implant




ProGuide™
Tri-Ball Tuner

FlowGuard™
Valved Introducer

ProGuide™
Long-Term Hemodialysis Catheter

> choose performance

Merit's **STEP-TIP** catheter has a large internal lumen and apposition deflection feature to minimize recirculation and optimize flow. This catheter addresses the needs of clinicians using the **OVER-THE-WIRE** technique and those who prefer the valved, **PEEL AWAY SHEATH** insertion technique.



CentrosFLO®
Long-Term Hemodialysis Catheter

FlowGuard™
Valved Introducer

CentrosFLO® Tunneler

implant

> improve patency

CentrosFLO Long-Term Hemodialysis Catheters are the next generation, **SELF-CENTERING** split-tip dialysis catheters designed by Dr. Stephen Ash. The CentrosFLO's self-centering design was created to address the issue of fibrin sheath formation. The catheter has **SMOOTH** surfaces and two preshaped, **CURVED TIPS** to automatically "center" the catheter within the vein.

Slip-Not®
Suture Retention Device



basixTOUCH™
Inflation Device



Fountain®
Infusion System




Squirt®
Fluid Dispensing System



treat



Prelude®
Short Sheath Introducer



ONE Snare®
Endovascular Snare System

EN Snare®
Endovascular Snare System

➤ treat effectively

The Prelude Short Sheath provides single access to dialysis fistulas and grafts for **DECLOT**, angioplasty, thrombolysis and temporary hemodialysis. Two snare configurations give you the options you need to strip fibrin sheath formation. A complete line of therapeutic **INFUSION** systems allows you to achieve more uniform dispersion of therapeutic agents. Merit's latest innovation in inflation device technology offers increased volume and pressure capability for **INFLATION**. The Slip-Not **SUTURE RETENTION** device helps manage bleeding after an AV fistula intervention by securing purse-string sutures.

treat

integrate





➤ customize for efficiency

To help you optimize the economics of providing dialysis quickly, efficiently, and safely, we offer **INTEGRATED** dialysis procedure packs with all of the products needed for your procedure. Having all of the products you need packaged together enables you to offer both **CONVENIENT** and **COST-EFFECTIVE** care.

peritoneal dialysis

prepare

place


select

integrate

assist

complete





prepare



Implantation Stencils

➤ plan for options

The **IMPLANTATION** stencil is designed to improve catheter placement **ACCURACY** by helping you select the optimal catheter style, size, implantation site and exit site for each patient. The stencil is conveniently included in every catheter kit and can be used for both left and right-sided placements. For upper abdomen and upper chest placements, Merit offers a complete set of three sterile **STENCILS** in each ExXTended catheter kit.



Y-TEC® VP-211
Implantation System

place


➤ implant with simplicity

Merit offers a comprehensive implantation **KIT** for the single-site incision of a PD catheter, avoiding the need to place the catheter through an over-sized laparoscopic port. This complete set of components is designed to streamline the implantation procedure and promote **EFFICIENCY**. The kit includes a unique cuff implantor, dilators, trocar assembly, and a tunneler. Multiple kit **OPTIONS** are available to support your specific implantation approach and tunneling style.

select




Flex-Neck® ARC™
Peritoneal Dialysis Catheter



Flex-Neck Classic
Peritoneal Dialysis Catheter





Flex-Neck ExxTended™
Peritoneal Dialysis Catheter

select

> choose with confidence

The Merit adult PD catheter has an internal diameter of 3.5 mm, allowing **35% HIGHER FLOW** rates than other catheters currently on the market. Catheters are available in various **SHAPES**, and tip-to-cuff **LENGTHS** to accommodate most patient types including infant, pediatric, adolescent, and small, medium, and large adults. Our Classic and ARC™ shapes are designed to accommodate various tunneling techniques and exit site options. Our complete ExxTended catheter kit includes all of the necessary tools for upper abdomen and upper chest (pre-sternal) procedures.



integrate

Gentle Lubricating Gel is a trademark of Pharmaceutical Innovations, Inc.
CUT HERE




integrate

➤ customize for efficiency


Our **COMPLETE** ExxTended catheter kit includes all of the necessary tools for upper abdomen and upper chest (pre-sternal) procedures. The kit includes a measuring tape, a tunneling tool, three stencils for precise patient marking, and a measuring rod to gauge the rectus-to-skin distance for more **ACCURATE** catheter sizing and cutting.

A stainless steel surgical instrument with a textured handle and a sharp, pointed tip.

Faller Trocar
Stainless Steel

A long, thin, black surgical instrument with a textured handle and a sharp tip.

**Implantation
Stylette**

A collection of small, clear plastic components used for repairing a PD catheter, including a long tube and several connectors.

Flex-Neck
PD Repair Kit

➤ fit with precision

Merit's PD accessories are specifically designed to **SUPPORT** catheter implantation, connection, and repair. The **REPAIR** Kit is designed to extend the life of the catheter. The metal Faller Trocar creates a clean, tight exit site, utilizing an antegrade tunneling approach. Three lengths of reusable, stainless steel implantation stylettes are specifically sized to **FIT** the interior diameter of our PD catheters.



Embedding® Tool

**Catheter
Connectors**

complete

➤ finalize with certainty

The Embedding Tool is used for dialysis **PREPLANNING** in cases where residual kidney function exists or dialysis will be needed in the future. This tool allows the catheter to be embedded under the skin during the initial implantation procedure. The proximal end of the catheter can be explanted later for future use, once the patient is **READY** for dialysis. This allows the tissue ingrowth cuffs to adhere to the muscle and subcutaneous tissue **IN ADVANCE** of dialysis. Titanium connectors, plastic replacement connectors, and caps are available as needed by the PD clinic or when the catheter is placed.



HEMODIALYSIS

prepare
access
navigate
implant
treat
integrate

PERITONEAL DIALYSIS

prepare
place
select
integrate
assist
complete

efficient. effective. safe.

At every step in your dialysis procedure, Merit Medical offers high-quality products designed to meet your clinical needs. From preparation to intervention, our products will help you achieve the clinical, economical, and procedural efficiency that is essential for delivering the best possible care.



Merit Medical Systems, Inc.

1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com

**Merit Medical Europe, Middle East,
and Africa (EMEA)**

Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria	0800 295 374	Luxembourg	8002 25 22
Belgium	0800 72 906 (Dutch) 0800 73 172 (French)	Netherlands	0800 022 81 84
Denmark	80 88 00 24	Norway	800 11629
France	0800 91 60 30	Sweden	020 792 445
Finland	0800 770 586	UK 0800	973 115
Germany	0800 182 0871		
Ireland (Republic)	1800 553 163		
Italy	800 897 005		

