



### Arterial Graft Component


6mm (ID) x 50cm

### Connector

6mm - 5mm (ID)

### Venous Outflow Component

5mm (ID), 6.3mm (OD), 19F (OD) x 40cm (customizable length)



Accessory Component Kit

### HeRO Graft Components

- Venous Outflow Component (HERO 1001)
- Arterial Graft Component (HERO 1002)
- Accessory Component Kit (HERO 1003)

### Other Possible Surgical Instruments (may be required, but not limited to the following list)

- |   |   |
|---|---|
| <input type="checkbox"/> C-Arm (for x-ray)*   | <input type="checkbox"/> 4x4 sterile gauze pads   |
| <input type="checkbox"/> Vascular ultrasound*   | <input type="checkbox"/> Radiographic contrast fluid  |
| <input type="checkbox"/> 5F micro-puncture set and 18 gauge needle*                     | <input type="checkbox"/> Various subcutaneous tissue & skin sutures   |
| <input type="checkbox"/> Various 0.035" guidewires<br>(stiff and floppy at least 150cm) | <input type="checkbox"/> Heparinized saline   |
| <input type="checkbox"/> Tissue tunneler set with 6mm and 7mm bullet tips               | <input type="checkbox"/> 6mm x 4cm x 80cm angioplasty balloon*  |
| <input type="checkbox"/> Various atraumatic vascular clamps                             | <input type="checkbox"/> 8mm x 4cm x 80cm angioplasty balloon*  |
| <input type="checkbox"/> Heavy duty scissors  | <input type="checkbox"/> Indeflator (for angioplasty balloon)*  |
| <input type="checkbox"/> Standard vessel loops  | <input type="checkbox"/> Various 5F catheters*  |
| <input type="checkbox"/> Syringe and syringe adapter                                    | <input type="checkbox"/> Hemodialysis catheter (for use during the bridging<br>period while ePTFE portion of the HeRO Graft<br>incorporates)* |

\*Other items to be considered which are NOT included in the Instructions for Use.

Surgical technique is at the discretion of the surgeon. Variations in technique and practices will inevitably and appropriately occur when clinicians take into account the needs of the individual patients, available resources, and limitations unique to an institution or type of practice.

Learn more at [MeritEMEA.com/hero](http://MeritEMEA.com/hero)



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

Merit Medical Europe, Middle East,  
& Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733