



CARDIAC INTERVENTION

# Hydrophilic Prelude SNAP™

Splittable Sheath Introducer

Now with  
**Hydrophilic**  
Coating

Requires  
**50%**  
less insertion  
force



**Easier to Break**  
requiring 66% less  
break force than  
its predecessor



**Seamless**  
sheath-to-dilator  
transition with  
hydrophilic coating



# Hydrophilic Prelude SNAP™

Splittable Sheath Introducer

The Hydrophilic Prelude SNAP features a **HYDROPHILIC** sheath coating that reduces insertion force by 50%\*. A smooth sheath to dilator transition can help reduce vessel trauma and improve smooth navigation through torturous anatomy. The Hydrophilic Prelude SNAP contains the same **HIGH QUALITY** hub which requires 66% less break force than its predecessor making it **EFFORTLESS** to remove and helping ensure leads stay in place.

## PRODUCT FEATURES

### Robust Valve Design

Tight seal to protect against air embolism and blood loss.

### Improved Grip

Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.




### Easier to Break

Requires 66% less break force than its predecessor.\*

### Seamless Transition

Smoother sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.



Competitors' sheath-to-dilator transitions

Our seamless sheath-to-dilator transition

### Hydrophilic Coating

Covers 75% of distal sheath length for ease of insertion while maintaining stability.

## KIT CONTENTS

- Splittable Sheath
- Dilator
- Guide Wire 0.038" x 50 cm or 0.038" x 80 cm
- 18G Needle
- 12 mL Syringe

## ORDERING INFORMATION

| Catalog Number with Sideport | French Size | Sheath Length | Hub Color     |
|------------------------------|-------------|---------------|---------------|
| PLSH-1006                    | 6           | 13 cm         | Green         |
| <b>PLSH-1007</b>             | <b>7</b>    | <b>13 cm</b>  | <b>Orange</b> |
| PLSH-1008                    | 8           | 13 cm         | Blue          |
| PLSH-1008.5                  | 8.5         | 13 cm         | Light Blue    |
| <b>PLSH-1009</b>             | <b>9</b>    | <b>13 cm</b>  | <b>Black</b>  |
| PLSH-1009.5                  | 9.5         | 13 cm         | Gray          |
| PLSH-1010                    | 10          | 13 cm         | Fuchsia       |
| PLSH-1010.5                  | 10.5        | 13 cm         | Light Fuchsia |
| PLSH-1011                    | 11          | 13 cm         | Yellow        |
| PLSH-1012                    | 12          | 13 cm         | Brown         |
| PLSH-1012.5                  | 12.5        | 13 cm         | Light Brown   |
| PLSH-2506                    | 6           | 25 cm         | Green         |
| PLSH-2507                    | 7           | 25 cm         | Orange        |
| PLSH-2508                    | 8           | 25 cm         | Blue          |
| PLSH-2509                    | 9           | 25 cm         | Black         |
| PLSH-2510                    | 10          | 25 cm         | Fuchsia       |
| PLSH-2510.5                  | 10.5        | 25 cm         | Light Fuchsia |
| PLSH-2511                    | 11          | 25 cm         | Yellow        |

| Catalog Number without Sideport | French Size | Sheath Length | Hub Color     |
|---------------------------------|-------------|---------------|---------------|
| PLSHX-1006                      | 6           | 13 cm         | Green         |
| PLSHX-1007                      | 7           | 13 cm         | Orange        |
| PLSHX-1008                      | 8           | 13 cm         | Blue          |
| PLSHX-1008.5                    | 8.5         | 13 cm         | Light Blue    |
| PLSHX-1009                      | 9           | 13 cm         | Black         |
| PLSHX-1009.5                    | 9.5         | 13 cm         | Gray          |
| PLSHX-1010                      | 10          | 13 cm         | Fuchsia       |
| PLSHX-1010.5                    | 10.5        | 13 cm         | Light Fuchsia |
| PLSHX-1011                      | 11          | 13 cm         | Yellow        |
| PLSHX-1012                      | 12          | 13 cm         | Brown         |
| PLSHX-1012.5                    | 12.5        | 13 cm         | Light Brown   |
| PLSHX-2506                      | 6           | 25 cm         | Green         |
| PLSHX-2507                      | 7           | 25 cm         | Orange        |
| PLSHX-2508                      | 8           | 25 cm         | Blue          |
| PLSHX-2509                      | 9           | 25 cm         | Black         |
| PLSHX-2510                      | 10          | 25 cm         | Fuchsia       |
| PLSHX-2510.5                    | 10.5        | 25 cm         | Light Fuchsia |
| PLSHX-2511                      | 11          | 25 cm         | Yellow        |

\*Data on file.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Merit.com

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 225180252  
(Français)  
+41 225948000  
(Italiano)  
+41 225180035  
**UK**  
0800 973 115