


# INTERVENTIONAL CRT SYSTEM EQUIPMENT LIST

Advancing LV lead delivery by simplifying straightforward cases and providing options for challenging anatomies



POWERED BY



[INTERVENTIONALCRT.COM](http://INTERVENTIONALCRT.COM)

## INTERVENTIONAL CRT SYSTEM EQUIPMENT LIST

### ACCESS

| DESCRIPTION                       | MANUFACTURER  | ORDER NUMBER                                                                                            | USES                                                                                                                                                                      | QUANTITY |
|-----------------------------------|---------------|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>PRELUDE SNAP™</b>              |               |                                                                                                         |                                                                                                                                                                           |          |
| 13 cm, with side port             | Merit Medical | 6F: PLS-1006<br>7F: PLS-1007<br>9F: PLS-1009<br>Available in 6F-12.5F<br>Available without side port    | • Introduction of RA or RV leads                                                                                                                                          |          |
| <b>HYDROPHILIC PRELUDE SNAP™</b>  |               |                                                                                                         |                                                                                                                                                                           |          |
| 13 cm, with side port             | Merit Medical | 6F: PLSH-1006<br>7F: PLSH-1007<br>9F: PLSH-1009<br>Available in 6F-12.5F<br>Available without side port | • Introduction of RA or RV leads when reduced insertion force is desired                                                                                                  |          |
| 25 cm, with side port             | Merit Medical | 6F: PLSH-2506<br>7F: PLSH-2507<br>9F: PLSH-2509<br>Available in 6F-11F<br>Available without side port   | • Introduction of RA or RV leads when reduced insertion force is desired                                                                                                  |          |
| <b>STIFFENED MINI ACCESS KIT™</b> |               |                                                                                                         |                                                                                                                                                                           |          |
| 5F, 15 cm                         | Merit Medical | S-MAK501N-15                                                                                            | • Especially useful during access when subclavian venoplasty is necessary<br>• Stiffened dilator provides additional support through scar tissue and/or difficult anatomy |          |

### LV LEAD DELIVERY SYSTEM

| DESCRIPTION                                           | MANUFACTURER  | ORDER NUMBER                                           | USES                                                                                                                                                       | QUANTITY |
|-------------------------------------------------------|---------------|--------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>WORLEY™ ADVANCED CORONARY SINUS GUIDE (CSG)</b>    |               |                                                        | <b>CS CANNULATION</b><br>• Provides 9F CS access and ability to perform Interventional CRT techniques                                                      |          |
| Standard, 40cm                                        | Merit Medical | WOR-CSG-B1-09                                          | • This size is used in majority of cases                                                                                                                   |          |
| Jumbo, 50cm                                           | Merit Medical | WOR-CSG-B2-09                                          | • Large curve for marked atrial enlargement                                                                                                                |          |
| Long, 50cm                                            | Merit Medical | WOR-CSG-BL1-09                                         | • May be useful with cephalic cut down or when deep seating and end support is needed                                                                      |          |
| Right Sided, 40cm                                     | Merit Medical | WOR-CSG-BR1-09                                         | • For CS access from the right side                                                                                                                        |          |
| <b>WORLEY™ ADVANCED LATERAL VEIN INTRODUCER (LVI)</b> |               |                                                        | <b>NAVIGATION &amp; SUB-SELECTION</b><br>• 7F allows for direct delivery of available LV leads<br>• 5.5F allows for direct delivery of LV leads up to 5.0F |          |
| Renal                                                 | Merit Medical | 7F: WORLVI-75-5-62-07-RE<br>5.5F: WORLVI-75-5-62-55-RE | • Renal Shape is used in majority of cases                                                                                                                 |          |
| Hockey Stick                                          | Merit Medical | 7F: WORLVI-75-5-62-07-HS<br>5.5F: WORLVI-75-5-62-55-HS | • Alternate shape for varying anatomies                                                                                                                    |          |
| Hook                                                  | Merit Medical | 7F: WORLVI-75-5-62-07-HO<br>5.5F: WORLVI-75-5-62-55-HO | • Alternate shape for varying anatomies                                                                                                                    |          |
| Multi-Purpose                                         | Merit Medical | 7F: WORLVI-75-5-62-07-MP<br>5.5F: WORLVI-75-5-62-55-MP | • Alternate shape for varying anatomies                                                                                                                    |          |
| <b>WORLEY™ VEIN SELECTOR</b>                          |               |                                                        | <b>NAVIGATION &amp; SUB-SELECTION</b>                                                                                                                      |          |
| Standard                                              | Merit Medical | 57538CS-WOR                                            | • Standard shape is used in majority of cases                                                                                                              |          |
| Vertebral                                             | Merit Medical | 57538CSV-WOR                                           | • Particularly useful with difficult CS access, navigation through Valve of Vieussens or Middle Cardiac Vein access                                        |          |
| Hook                                                  | Merit Medical | 57538CSHK-WOR                                          | • Useful with extreme acute target branch take off angles                                                                                                  |          |
| <b>WORLEY™ ADVANCED CONTRAST ADMINISTRATION KIT</b>   |               |                                                        | <b>VISUALIZATION</b>                                                                                                                                       |          |
| With Contrast Bowl                                    | Merit Medical | K12-WORLEY1                                            | • For easy administration of contrast by assistant while allowing wire use                                                                                 |          |
| Without Contrast Bowl                                 | Merit Medical | K12-WORLEY2                                            | • For easy administration of contrast by assistant while allowing wire use                                                                                 |          |

## INTERVENTIONAL CRT SYSTEM EQUIPMENT LIST

### VENOPLASTY

| DESCRIPTION                                      | MANUFACTURER      | ORDER NUMBER                                               | USES                                                                                                                                                | QUANTITY |
|--------------------------------------------------|-------------------|------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>0.035" LUMEN SUBCLAVIAN BALLOON EXAMPLES*</b> |                   |                                                            | <b>*Many manufacturer options are available. Key characteristics include: Ultra non-compliant, over the wire, 0.035" lumen, RBP at least 14 ATM</b> |          |
| Conquest®                                        | Bard              | 9mm x 4cm: CQ7594<br>6mm x 4cm: CQ7564                     | • For subclavian venoplasty when 0.035" wire is across occlusion                                                                                    |          |
| Mustang® Balloon Dilatation Catheter             | Boston Scientific | 9mm x 4cm: 3917109047<br>6mm x 4cm: 3917106047             | • For subclavian venoplasty when 0.035" wire is across occlusion                                                                                    |          |
| <b>0.018" LUMEN SUBCLAVIAN BALLOON EXAMPLE*</b>  |                   |                                                            | <b>*Many manufacturer options are available. Key characteristics include: Ultra non-compliant, over the wire, 0.018" lumen, RBP at least 14 ATM</b> |          |
| NC Sterling®                                     | Boston Scientific | 6mm x 4 cm:<br>UPN: H74939032604040<br>GTIN: 8714729123781 | • For subclavian venoplasty when 0.018" wire is across occlusion                                                                                    |          |
| <b>IMPRESS® HYDROPHILIC CATHETER</b>             |                   |                                                            |                                                                                                                                                     |          |
| 5F, 65cm, 0.046" ID, Angled Tip                  | Merit Medical     | 56538KA2-H                                                 | • To help direct wire across subclavian occlusion or to upgrade to larger wire                                                                      |          |
| 4F, 65cm, 0.040" ID, Straight Tip                | Merit Medical     | 46538STS-H                                                 | • To exchange 0.035" glide wire for 0.035" extra stiff support wire                                                                                 |          |
| <b>CORONARY BALLOON EXAMPLES*</b>                |                   |                                                            | <b>*Many manufacturer options are available. Key characteristics include: Non-compliant, rapid exchange, RBP 18-20 ATM</b>                          |          |
| NC Sprinter®                                     | Medtronic         | 3mm x 15mm: NCSP3015X                                      | • Coronary vein venoplasty                                                                                                                          |          |
| NC Trek®                                         | Abbott            | 3mm x 15mm: 1012449-15                                     | • Coronary vein venoplasty                                                                                                                          |          |
| <b>basixTOUCH™ INFLATION DEVICE</b>              |                   |                                                            | <b>VISUALIZATION</b>                                                                                                                                |          |
| 35 ATM/BAR                                       | Merit Medical     | IN8130                                                     | • Inflation of subclavian or coronary balloons                                                                                                      |          |

### LV LEAD PLACEMENT WITH SNARING

| DESCRIPTION                              | MANUFACTURER  | ORDER NUMBER | USES                                                                                                                                                                | QUANTITY |
|------------------------------------------|---------------|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>ONE SNARE®</b>                        |               |              |                                                                                                                                                                     |          |
| 10mm loop, 120cm, 4F                     | Merit Medical | ONE1000      | • Snaring of LV lead delivery wire to create rail for lead placement in small or tortuous vessels                                                                   |          |
| <b>MERIT SURECROSS® SUPPORT CATHETER</b> |               |              |                                                                                                                                                                     |          |
| 1.52F crossing profile, 90cm             | Merit Medical | SC1490       | • To protect collateral vessels when pulling 0.014" wire through<br>• Can be useful to help navigate collateral vessels when more support and pushability is needed |          |
| <b>IMPRESS® HYDROPHILIC CATHETER</b>     |               |              |                                                                                                                                                                     |          |
| 4F, 65cm, 0.040" ID, KA2 Tip             | Merit Medical | 46538KA2-H   | • To direct 0.035" wire into exit vein for subsequent snare catheter placement over the wire. To then be replaced by snare                                          |          |

CUSTOMER SERVICE 1 800 35 MERIT **MERIT.COM**



# INTERVENTIONAL CRT SYSTEM EQUIPMENT LIST

## GUIDE WIRES

| DESCRIPTION                                                | MANUFACTURER      | ORDER NUMBER                                     | USES                                                                                                                           | QUANTITY |
|------------------------------------------------------------|-------------------|--------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>0.035" INQWIRE™ DIAGNOSTIC GUIDE WIRE</b>               |                   |                                                  |                                                                                                                                |          |
| Standard, 180cm, Jtip                                      | Merit Medical     | IQ35F180J3                                       | • Provides support for CSG in CS                                                                                               |          |
| Firm, 180cm, Jtip                                          | Merit Medical     | IQH35F180J3F                                     | • Provides support for CSG in CS                                                                                               |          |
| <b>MERIT LAUREATE® HYDROPHILIC GUIDE WIRE</b>              |                   |                                                  |                                                                                                                                |          |
| 0.018", 180cm, Angled Tip, Polymer Coated                  | Merit Medical     | LWSTDA18180                                      | • Crossing subclavian occlusions<br>• Venoplasty balloon delivery                                                              |          |
| 0.035", 180cm, Angled Tip, Polymer Coated                  | Merit Medical     | LWSTDA35180                                      | • Crossing subclavian occlusions<br>• Venoplasty balloon delivery<br>• Facilitating CS cannulation                             |          |
| <b>LEAD PLACEMENT WIRE EXAMPLES*</b>                       |                   |                                                  | <b>*Many manufacturer options are available. Key characteristics include: Polymer Tip (PT), Straight Tip, Length, Diameter</b> |          |
| CholCE PT™ Floppy<br>0.014", Straight Tip                  | Boston Scientific | 182cm: H749 12160-01 2<br>300cm: H749 12154-01 2 | • Navigation of collaterals<br>• Lead delivery                                                                                 |          |
| CholCE PT™ Extra Support<br>0.014", Straight Tip           | Boston Scientific | 182cm: H749 12161-01 2<br>300cm: H749 12155-01 2 | • Navigation of collaterals<br>• Buddy wire support<br>• Lead delivery<br>• Antidromic Snare Technique                         |          |
| V-18™ ControlWire™<br>0.018", Polymer Sleeve, Straight Tip | Boston Scientific | 200cm: M001468520                                | • Buddy wire/rail support<br>• Medtronic LV lead delivery                                                                      |          |
| <b>TORQUE DEVICE</b>                                       |                   |                                                  | <b>TORQUING OF WIRES</b>                                                                                                       |          |
| Sea Dragon™                                                | Merit Medical     | TD500                                            | • Fits 0.018"-0.038" guide wires                                                                                               |          |
| Merit                                                      | Merit Medical     | MAP500                                           | • Fits 0.009"-0.018" guide wires                                                                                               |          |

