

An astronaut in a full space suit is shown in a spacewalk, floating against a backdrop of a starry galaxy and the curved horizon of Earth. The astronaut's helmet visor is highly reflective, showing a bright light source. The overall scene is dramatic and emphasizes reaching new frontiers.


**REACH
BEYOND
LIMITS**



Merit Pursue™

Microcatheter

For the distal vasculature
that seemed out-of-reach.



Catheter Compatibility

French Size	Micro-catheter ID	Minimum Guide Catheter ID	Maximum Guide Wire OD	Maximum Particle Size	Maximum Micro-sphere Size	Maximum Coil Size
2.8F – 1.7F	0.016" (0.40 mm)	0.040" (1.02 mm)	0.014" (0.36 mm)	≤ 500 μm	≤ 500 μm	0.014" (0.36 mm)
2.9F – 2.0F	0.020" (0.50 mm)	0.042" (1.07 mm)	0.018" (0.46 mm)	≤ 710 μm	≤ 700 μm	0.018" (0.46 mm)

The Merit Pursue™ Microcatheter is a small distal French size microcatheter designed for pushability and trackability through small and tortuous vessels.

The small outer diameter allows you to get deeper into the vasculature and still deliver the same spherical and non-spherical embolics of larger diameter microcatheters.

Chemical Compatibility

- Ethanol
- Lipiodol
- DMSO (Dimethyl Sulfoxide)
- Cisplatin
- Doxorubicin
- Irinotecan
- Cyanoacrylate

Ordering Information

Catalog Number	Proximal OD	Distal OD	Usable Length (cm)	Tip Shape	Hydrophilic Coating Length (cm)	Distance From Tip to Markerband (mm)	Priming Vol. (mL)
28HC1711045	2.8F	1.7F	110	45 degree	80	0.6	0.42
28HC17110SN	2.8F	1.7F	110	Swan	80	0.6	0.42
28HC17110ST	2.8F	1.7F	110	Straight	80	0.6	0.42
28HC1713045	2.8F	1.7F	130	45 degree	80	0.6	0.50
28HC17130SN	2.8F	1.7F	130	Swan	80	0.6	0.50
28HC17130ST	2.8F	1.7F	130	Straight	80	0.6	0.50
28HC1715045	2.8F	1.7F	150	45 degree	80	0.6	0.53
28HC17150SN	2.8F	1.7F	150	Swan	80	0.6	0.53
28HC17150ST	2.8F	1.7F	150	Straight	80	0.6	0.53
29HC2011045	2.9F	2.0F	110	45 degree	80	0.6	0.50
29HC20110SN	2.9F	2.0F	110	Swan	80	0.6	0.50
29HC20110ST	2.9F	2.0F	110	Straight	80	0.6	0.50
29HC2013045	2.9F	2.0F	130	45 degree	80	0.6	0.57
29HC20130SN	2.9F	2.0F	130	Swan	80	0.6	0.57
29HC20130ST	2.9F	2.0F	130	Straight	80	0.6	0.57
29HC2015045	2.9F	2.0F	150	45 degree	80	0.6	0.63
29HC20150SN	2.9F	2.0F	150	Swan	80	0.6	0.63
29HC20150ST	2.9F	2.0F	150	Straight	80	0.6	0.63

Contact your Merit Medical sales representative for more information or visit our website www.merit.com/pursue.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 91 1238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115