




# PATIENT SAFETY IS NOT NEGOTIABLE

**DualCap®**  
Disinfection & Protection System

**DTXPlus®  
& Safedraw®**  
Disposable Transducer and  
Closed in-line Sampling System



# REDUCING MEDICAL ERRORS

**DualCap® and DTXPlus® with Safedraw® products – engineered with several innovative features designed to prevent medical errors.**

## **Straight Fluid Column**


Clear housing easily identifies air in the line which can lead to damping and inaccurate information that guides decision making.

## **Safedraw Sample Port**

No dead space, minimizes air, damping and inaccurate information that can result in errors.

## **The DualCap System Locks In Protection**

The only cap for male luer connectors that completely blocks disinfectant from the patient's bloodstream<sup>1</sup>.



## **Safedraw Syringe**

Closed draw system to reduce line contamination and unwanted air, critical features designed to reduce medical errors.

## **EasyVent® Luer Lock**

Zero and prime without touch contamination assisting in the effort to protect against infection.

## **Angled Cannula**

Minimizes air, reduces damping designed for improved accuracy to assist clinicians taking action.

Catalog Number	Description	Length	Tubing Color	Units Per Box/Carton
SAFEDRAW® CLOSED IN-LINE SAMPLING SYSTEM				
682142	Safedraw® with Meritans® DTXPlus® 60" Pressure Monitoring Kit and Blood Sampling System	120 / 30 cm	Clear	5/25
682143	Safedraw® with Meritans® DTXPlus® 72" Pressure Monitoring Kit and Blood Sampling System	150 / 30 cm	Clear	5/25

## SAFEDRAW® STAND-ALONE SETS

682144	Safedraw® 60" Blood Sampling System	120/30 cm	Clear	5/50
682145	Safedraw® 72" Blood Sampling System	150/30 cm	Clear	5/50

## SAFEDRAW® ACCESSORIES

682158	Safedraw® for BD Vacutainer® Transfer Device	N/A	N/A	25/125
682159	Safedraw® for BD Luer-Lok™ Safe-Needle	N/A	N/A	50/500

Catalog Number	Product Description	Caps Per Box	Boxes Per Case	Total Caps Per Case
DUALCAP® DISINFECTION & PROTECTION SYSTEM				
001-300	DualCap Solo™ Cap with foil lid for needle-free valve, packed in dispenser box.	350	12	4200
001-400	DualCap Solo™ Cap with foil lid for male luer connector, packed in dispenser box.	200	12	2400
500-LB	DualCap® IV Pole Strips® For needle-free valves. 35 strips/10 caps per strip.	350	12	4200
600-DB	DualCap® IV Pole Strips® For male luer connectors. 35 strips/5 caps per strip.	175	6	1050

<sup>1</sup> Avomeen Analytical Services, Testing Report prepared for Catheter Connections, October 2013.

**Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.**



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

[merit.com](http://merit.com)

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115