




KNOW YOUR DIRECTION



HeartSpan®
Steerable Sheath Introducer



NEUTRAL POSITION




INDICATOR

HeartSpan®

Steerable Sheath Introducer

HIGHER CLINICAL SUCCESS RATE*



The study showed that the use of a steerable sheath improved catheter stability and catheter-to-tissue contact, which resulted in higher clinical treatment success for AF ablation.*

KIT CONTENTS

- (1) 8.5 F Steerable Sheath Introducer
- (1) Dilator
- (1) 0.032" x 180 cm guidewire

ORDERING INFORMATION

Part Number	Curl Type	Curl Size	Compatible Needle Length
HST-085-00	Small	16.4 mm	98 cm
HST-085-10	Medium	22.4 mm	98 cm
HST-085-20	Large	50.0 mm	98 cm

*Piorkowski, Christopher, et al. "Steerable Versus Nonsteerable Sheath Technology in Atrial Fibrillation Ablation, A Prospective, Randomized Study" *Circulation: Arrhythmia and Electrophysiology* 4.2 (2011): 157-165.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
 1600 West Merit Parkway
 South Jordan, Utah 84095
 1.801.253.1600
 1.800.35.MERIT

merit.com

Merit Medical Europe, Middle East, & Africa (EMEA)
 Amerikalaan 42, 6199 AE Maastricht-Airport
 The Netherlands
 +31 43 358 82 22

Merit Medical Ireland Ltd.
 Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733

Australia
 +61 (0) 3 8373 4943

Austria
 0800 295 374

Belgium
 0800 72 906 (Dutch)
 0800 73 172 (Français)

Brazil
 +55 11 4561 - 0062

Canada
 1.800.364.4370

China
 +86-10 8561 0788

Denmark
 80 88 00 24

Finland
 0800 770 586

France
 0800 91 60 30

Germany
 0800 182 0871

Hong Kong
 +852-2207 6888

India
 +91-80-41223376

Ireland (Republic)
 1800 553 163

Italy
 800 897 005

Luxembourg
 8002 25 22

Netherlands
 0800 022 81 84

Norway
 800 11629

Portugal
 308 801 034

Russia
 +7 495 221 89 02

Singapore
 +65 67536311

South Korea
 +8210 2205 0707

Spain
 +34 911238406

Sweden
 020 792 445

Switzerland
 (Deutsch)
 +41 22 518 02 30
 (Français)
 +41 22 518 02 52


(Italiano)
 +41 22 518 00 35

UK
 0800 973 115

NEUTRAL POSITION INDICATOR


Know your Direction

With a tactile click at neutral position, move from 0° to 180° to achieve the desired curve within the atrial chamber.





LOWER CROSSING PROFILE

Seamless dilator to sheath transition for smooth advancement across the atrial septum.



ADDITIONAL FEATURES



Atraumatic tip




Braided construction



Distal side holes to minimize cavitation



Side port with three-way stopcock



Radiopaque marker band



Hemostatic valve that reduces bleed back and air embolism



Lubricious coating