

OPTIONS MATTER

Hemostasis Valves
& Accessories



Hemostasis Valves & Accessories



FLO40XR™

The FLO40XR hemostasis valve offers a shorter body and extension tubing as an additional feature. Valve rated up to 400 psi.

SPECIFICATIONS

Inner lumen: 0.096" (2.44mm)/7.3F
Extension tube length 15cm (6")



FLO50™

The FLO50 hemostasis valve offers a shorter body maintaining a large bore inner diameter rated up to 400 psi.

SPECIFICATIONS

Inner lumen: .118" (3.0mm)/9F



Honor®

The easy push-pull mechanism of the Honor hemostasis valve allows single-handed manipulation to open and close the valve while introducing devices into the vasculature.

SPECIFICATIONS

Inner lumen: .094" (2.4mm)/7.2F



DoublePlay™

Two large bore hemostasis ports allow easier manipulation and control of interventional devices during "kissing balloon" techniques.

SPECIFICATIONS

Main inner lumen: .120" (3.1mm)/9F



AccessPLUS™

The AccessPLUS large bore hemostasis valve has optimal sealing surface for blood control.

SPECIFICATIONS

Inner lumen: .120" (3.1mm)/9F



Passage®

The Passage represents the first of Merit's standard line of hemostasis valves.

SPECIFICATIONS

Inner lumen: .097" (2.5mm)/7.4F

Access the right hemostasis options



FLO30™

Merit Medical's FLO30 hemostasis valve is a short-body adapter providing maximum working length when required.

SPECIFICATIONS

Inner lumen: .118" (3mm)/9F



MBA™ - Merit's Bleedback Alternative™

A premier hemostasis valve which features two seals. A conventional "touhy" silicone hemostasis valve provides security that the interventional device will stay in place. The secondary slit seal delivers a virtually bloodless procedure.

SPECIFICATIONS

Inner lumen: .118" (3mm)/9F

Merit Accessory Pak

Merit provides multiple hemostasis valve accessory combinations tailored to your interventional needs.



Hemostasis valves and accessories are also available with any Merit inflation device. Kits are available in standard catalog numbers as well as custom kit configurations.

BasixTOUCH™ (pictured)
Blue Diamond™

BasixCOMPAK™
BasixTOUCH40™
DiamondTOUCH™



Ordering Information

Catalog Number	Hemostasis valve	Units per box
Standalone Hemostasis Valves		
MAP150	Access plus, large bore	25
MAP220	Double Play, large bore	25
MAP300	Honor, small bore	25
MAP350	Honor, small bore with bonded 20cm tubing and a stopcock	25
MAP400	MBA, large bore	25
MAP450	MBA, large bore with bonded 25.4cm high pressure tubing	25
MAP900	Standalone FLO50 Valve	25
FLO30	FLO30, large bore	25
FLO40XR	Standalone FLO40XR Valve	20
K12-00040	Passage, small bore with bonded 15.2cm high pressure tubing	25

Body material: Polycarbonate



Catalog Number	Hemostasis valve	Units per box	Accessories		
			Metal wire introducer	Plastic wire introducer	Torque device
Hemostasis Valve MAP Kits					
MAP101	Passage, small bore	25	X		
MAP102	Passage, small bore	25	X		X
MAP111	Access 9, large bore	25	X		
MAP112	Access 9, large bore	25	X		X
MAP113	Access 9, large bore	25		X	
MAP151	Access plus, large bore	25	X		
MAP152	Access plus, large bore	25	X		X
MAP153	Access plus, large bore	25		X	
MAP154	Access plus, large bore	25		X	X
MAP221	Double Play, large bore	25	X		
MAP222	Double Play, large bore	25	X		X
MAP301	Honor, small bore	25	X		
MAP302	Honor, small bore	25	X		X
MAP351	Honor, small bore with bonded 20cm tubing and a stopcock	25	X		
MAP352	Honor, small bore with bonded 20cm tubing and a stopcock	25	X		X
MAP353	Honor, small bore with bonded 50cm tubing and a stopcock	25	X		X
MAP403	MBA, large bore	25	X		X
MAP453	MBA, with bonded 25.4cm high pressure tubing	25	X		X
MAP952	FLO50 Valve, Metal Guide Wire Insertion Tool, & Torque Device	25	X		X
K05	Custom Kits				

Accessories

MAP500	Torque Device 0.018", salmon	25			X
MAP501	Torque Device 0.018", green	25			X
MAP550	Metal Wire Insertion Tool	25	X		
MAP600	Plastic Wire Insertion Tool	25		X	

Body material: Polycarbonate

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115