




SPINR™

High-Performance Guidewire Controller


Precision. Versatility. Adaptability.™



PRECISION CONTROL
Ranging from 0-2500 RPM



ADAPTABLE
Works with guidewires from
0.014" to 0.038"



VERSATILE
Approved for both coronary
& peripheral procedures

SPINR is used to maneuver guidewires in the coronary and peripheral vasculature during interventional or diagnostic procedures. SPINR is not intended for use in the neurovasculature.

SPINR™

High-Performance Guidewire Controller

The **SPINR** high-performance guidewire controller can help when accessing distal anatomy, navigating tortuous vessels, crossing old lesions, and working through fresh blockages.

Providing a wide range of versatility, the **SPINR** puts the control in your hands, so you, the clinician can easily and precisely adjust rotational speed.

- **Precision Control** – Manipulate guidewire at adjustable speeds from 0 to 2,500 RPM
- **Versatility** – Intended for a wide range of procedures in both coronary and peripheral vasculature
- **Adaptability** – Compatible with your preferred high-torque 0.014" - 0.038" guidewire
- **Time saver** – Faster than traditional torque devices
- **Cost saver** – Less expensive than comparable electromechanical devices



Getting You There Faster. Discover how the SPINR can help you get your guidewire where you want it.

Contact your **Merit Medical Sales Professional** to experience high-performance guidewire control with the **SPINR**.

Product Number	Description	Guidewire Compatibility
3X-GW	SPINR Controller 3 Clockwise & Counterclockwise Rotations per Actuation	0.014"-0.038"
5X-GW	SPINR Controller 5 Clockwise & Counterclockwise Rotations per Actuation	0.014"-0.038"



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30
(Français)
+41 22 518 02 52
(Italiano)
+41 22 518 00 35

UK
0800 973 115