



# SPEED MEETS PRECISION



## EmboCube™ Embolization Gelatin

EmboCube is the only gelatin foam that is uniformly cut and loaded into a syringe for faster preparation and delivery.



S P E E D • P R E C I S I O N • C O N S I S T E N C Y ™

# EMBOCUBE™ EMBOLIZATION GELATIN

## DELIVERY



- Green color coded EmboCube can be delivered through microcatheters with a  $\geq 0.024''$  (0.61 mm) inner diameter, including:



Maestro® Microcatheter,  
2.8F - 2.8F and 2.9F - 2.9F



EmboCath® Plus Microcatheter,  
3.0F - 2.9F



- Red color coded EmboCube can be delivered through catheters with a  $\geq 0.040''$  (1.02 mm) inner diameter.

## ORDERING INFORMATION

CATALOG NUMBER	HYDRATED SIZE	WEIGHT	RECOMMENDED INJECTION SYRINGE	MINIMUM CATHETER ID	COLOR CODE
EC2525	2.5 mm	25 mg	1 mL	0.024" (0.61 mm)	●
EC2550	2.5 mm	50 mg	1 mL	0.024" (0.61 mm)	●
EC25100	2.5 mm	100 mg	1 mL	0.024" (0.61 mm)	●
EC5025	5.0 mm	25 mg	3 mL	0.040" (1.02 mm)	●
EC5050	5.0 mm	50 mg	3 mL	0.040" (1.02 mm)	●
EC50100	5.0 mm	100 mg	3 mL	0.040" (1.02 mm)	●

Packaged 5/box.

## EMBOCUBE DELIVERY KIT

Each kit includes two injection syringes (1 mL and 3 mL) and a 3-way stopcock for added convenience.

<b>K12-09739A</b>	1 mL Medallion® Syringe	●	3 mL Medallion® Syringe	●	3-way Stopcock
-------------------	-------------------------	---	-------------------------	---	----------------

Packaged 10/box.

Contact your Merit Medical sales representative for more information or visit [www.merit.com/emboCube](http://www.merit.com/emboCube).

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

[merit.com](http://merit.com)

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115