





EXPANDING OPTIONS

FOR HEMODYNAMIC MONITORING



Merit Pulmonary Artery and
Thermodilution Catheters

Merit Pulmonary Artery and Thermodilution Catheters



Merit's Safety-Wedge

PULMONARY ARTERY CATHETERS

Product Code	French Size	Lumens	Length	Accepted Guide Wire	Units per Box
--------------	-------------	--------	--------	---------------------	---------------

PVC Material

497527	6F	2	110 cm	.025" Wire	5
497523	7F	2	110 cm	.025" Wire	5
497540	7F	2	110 cm	.035" Wire	5

Polyurethane Material

TD2502NXF	5F	2	90 cm	.018" Wire	5
TD2602NXF	6F	2	110 cm	.025" Wire	5
TD2702NXF	7F	2	110 cm	.025" Wire	5
CL2702NXF	7F	2	110 cm	.025" Wire	5
TD2703NF	7F	3	110 cm	.025" Wire	5

CRITIKIT™ ICE/ROOM TEMPERATURE INJECTATE SYSTEMS AND ACCESSORIES

Product Code	Description	Sterile	Box/Case Content (Units)
680000	Ice/Room temperature Injectate Critikit (with in-line sensor, control syringe & cooling coil)	Yes	4/8
680002	Injectate delivery system for performing cardiac output measurements.	Yes	5
680006	In-line injectate sensor housing	Yes	25/250
682252	10 mL control syringe	Yes	10/50

Monitoring Life™

- Utilize multiple lumens to meet advanced monitoring needs.
- Merit's unique Safety-Wedge prevents over-pressurization and potential pulmonary artery rupture.
- Contamination shield provides additional barrier for infection control.
- Available in 5–7.5 French catheters, with 2–5 lumens.

THERMODILUTION CATHETERS

Product Code	French Size	Lumens	Length	Accepted Guide Wire	Special Features	Units per Box
--------------	-------------	--------	--------	---------------------	------------------	---------------

PVC Material/Standard

497537	6F	4	110 cm	.025" Wire	Safety Wedge/Sleeve	5
497521	7F	4	110 cm	.025" Wire	Safety Wedge/Sleeve	5
497539	7F	4	110 cm	.025" Wire	Safety Wedge	5
497526	7F	4	110 cm	.025" Wire	-	5
497520	7.5F	5	110 cm	.025" Wire	Safety Wedge/Sleeve	5

PVC Material/S-Curve Tip

497534	7F	4	110 cm	.025" Wire	Safety Wedge/Sleeve	5
497536	7F	4	110 cm	.025" Wire	-	5

PVC Material/Soft Body

497531	7.5F	5	110 cm	.025" Wire	Safety Wedge/Sleeve	5
--------	------	---	--------	------------	---------------------	---

Polyurethane Material/Standard

TD2504NDXF	5F	4	90 cm	.018" Wire	Safety Wedge	5
TD2604NDXF	6F	4	110 cm	.025" Wire	Safety Wedge	5
TD2604NXF	6F	4	110 cm	.025" Wire	-	5
TD2704NDF	7F	4	110 cm	.025" Wire	Safety Wedge/Sleeve	5
TD2704NXF	7F	4	110 cm	.025" Wire	-	5
TD2755NDF	7.5F	5	110 cm	.025" Wire	Safety Wedge/Sleeve	5
TD2755NDPXF	7.5F	5	110 cm	.025" Wire	Safety Wedge/Paceport	5
TD2755NXF	7.5F	5	110 cm	.025" Wire	-	5

Polyurethane Material/S-Curve Tip

TD2704NCXF	7F	4	110 cm	.025" Wire	-	5
------------	----	---	--------	------------	---	---

Polyurethane Material/Soft Body

TD2755NX	7.5F	5	110 cm	.025" Wire	-	5
----------	------	---	--------	------------	---	---

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733