





EXPANDING OPTIONS

FOR HAEMODYNAMIC MONITORING



Merit Pulmonary Artery and
Thermodilution Catheters

Merit Pulmonary Artery and Thermodilution Catheters



Monitoring Life™

- Utilize multiple lumens to meet advanced monitoring needs.
- Merit's unique Safety-Wedge prevents over-pressurization and potential pulmonary artery rupture.
- Contamination shield provides additional barrier for infection control.
- Available in 5–7.5 French catheters, with 2–5 lumens.

THERMODILUTION CATHETERS

Product Code	French Size	Lumens	Length	Accepted Guide Wire	Special Features	Units per Box
Polyurethane Material/Standard						
TD2504NDF	5F	4	90 cm	.018" Wire	Safety Wedge	5
TD2604NDF	6F	4	110 cm	.025" Wire	Safety Wedge	5
TD2604NXF	6F	4	110 cm	.025" Wire	-	5
TD2704NDF	7F	4	110 cm	.025" Wire	Safety Wedge/Sleeve	5
TD2704NXF	7F	4	110 cm	.025" Wire	-	5
TD2755NDF	7.5F	5	110 cm	.025" Wire	Safety Wedge/Sleeve	5
TD2755NDPXF	7.5F	5	110 cm	.025" Wire	Safety Wedge/Paceport	5
TD2755NXF	7.5F	5	110 cm	.025" Wire	-	5
Polyurethane Material/S-Curve Tip						
TD2704NCF	7F	4	110 cm	.025" Wire	-	5
Polyurethane Material/Soft Body						
TD2755NX	7.5F	5	110 cm	.025" Wire	-	5

PULMONARY ARTERY CATHETERS

Product Code	French Size	Lumens	Length	Accepted Guide Wire	Units per Box
Polyurethane Material					
TD2502NXF	5F	2	90 cm	.018" Wire	5
TD2602NXF	6F	2	110 cm	.025" Wire	5
TD2702NXF	7F	2	110 cm	.025" Wire	5
CL2702NXF	7F	2	110 cm	.025" Wire	5
TD2703NF	7F	3	110 cm	.025" Wire	5

CRITIKIT™ ICE/ROOM TEMPERATURE INJECTATE SYSTEMS AND ACCESSORIES

Product Code	Description	Sterile	Box/Case Content (Units)
680000	Ice/Room temperature Injectate Critikit (with in-line sensor, control syringe & cooling coil)	Yes	4/8
680006	In-line injectate sensor housing	Yes	25/250
680178	Ice/Room temperature Injectate Critikit (with in-line sensor and control Syringe, without cooling coil)	Yes	5/20
682252	10 mL control syringe	Yes	10/50

CRITICATH® PULMONARY ARTERY CATHETERS

Product Code	Model	Size	Length	No. of Lumens	PA Wedge Pressure	CV Pressure	Cardiac Output	CVP Medic Port	Additional Information	Box/Case Content (units)
680372	SP5325	5F	80 cm	2	Yes	No	No	No	Monitoring Catheter	5/20
680373	SP5325L	5F	110 cm	2	Yes	No	No	No	Monitoring Catheter	5/20
680374	SP5327	7F	110 cm	2	Yes	No	No	No	Monitoring Catheter	5/20

CRITICATH® THERMODILUTION CATHETERS

Product Code	Model	Size	Length	No. of Lumens	PA Wedge Pressure	CV Pressure	Cardiac Output	CVP Medic Port	Additional Information	Box/Case Content (units)
680349	SP5105	5F	80 cm	4	Yes	Yes	Yes		PVC Catheter Material	5/20
680103	SP5107U	7F	110 cm	4	Yes	Yes	Yes			5/20
680381	SP5107U-14	7F	110 cm	4	Yes	No	No		Pre-mounted Luer-Lok™ AC shield	5/20
680382	SP5107U + SP5045	7F	110 cm	4	Yes	Yes	Yes		Pre-mounted In-Line Inject, sensor	5/20
680378	SP5507U	7.5F	110 cm	5	Yes	Yes	Yes	Yes		5/20
680379	SP5507U-14	7.5F	110 cm	5	Yes	Yes	Yes	Yes	Pre-mounted Luer-Lok™ AC shield	5/20
680380	SP5507U + SP5045	7.5F	110 cm	5	Yes	Yes	Yes	Yes	Pre-mounted In-Line Inject, sensor	

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733