


DuraFlow 2


CHRONIC DIALYSIS CATHETER



DuraFlow 2

CHRONIC DIALYSIS CATHETER

Engineered for Durability



Ordering Information (5 units per box)

Description	Basic Kit	VascPak Kit	Catheter Only
20 cm	<input type="checkbox"/> H787103012145	—	—
22 cm	<input type="checkbox"/> H787103012125	—	—
24 cm	<input type="checkbox"/> H787103012015	<input type="checkbox"/> H787103012185	<input type="checkbox"/> H787103017015
28 cm	<input type="checkbox"/> H787103012025	<input type="checkbox"/> H787103012195	<input type="checkbox"/> H787103017025
32 cm	<input type="checkbox"/> H787103012035	<input type="checkbox"/> H787103012205	<input type="checkbox"/> H787103017035
36 cm	<input type="checkbox"/> H787103012045	<input type="checkbox"/> H787103012215	<input type="checkbox"/> H787103017045
40 cm	—	—	<input type="checkbox"/> H787103017055
48 cm	—	—	<input type="checkbox"/> H787103017155
55 cm	—	—	<input type="checkbox"/> H787103017145

Ordering Information (1 unit per box)

Description	Basic Kit	VascPak Kit
24 cm	<input type="checkbox"/> H787103012011	—
28 cm	<input type="checkbox"/> H787103012021	—
32 cm	<input type="checkbox"/> H787103012031	—
36 cm	<input type="checkbox"/> H787103012041	—
40 cm	<input type="checkbox"/> H787103012055	<input type="checkbox"/> H787103012225
48 cm	<input type="checkbox"/> H787103012155	<input type="checkbox"/> H787103012235
55 cm	<input type="checkbox"/> H787103012115	<input type="checkbox"/> H787103012245

BASIC KIT INCLUDES: (1) DuraFlow 2 Dialysis Catheter, (1) 16F Valved Breakaway Introducer, (1) 12F Dilator, (1) 14F Dilator, (1) Tri-Ball Tunneler with Sleeve, (2) Injection Caps, (1) #11 Blade Safety Scalpel, (1) 0.038" J/Flex Guidewire, (1) Adhesive Wound Dressing, (1) 18 Ga. x 2-3/4" Introducer Needle

VASCPAK[®] KIT INCLUDES: (1) DuraFlow 2 Dialysis Catheter, (1) Adhesive Wound Dressing, (1) 16F Valved Peelable Sheath, (1) Tri-ball Tunneler with Sleeve, (2) Injection Caps, (1) 12F Vessel Dilator, (1) 14F Vessel Dilator

IMPORTANT RISK INFORMATION

INDICATIONS FOR USE: The DuraFlow 2 Hemodialysis Catheter is indicated for use in attaining long-term vascular access for hemodialysis and apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required. Catheters greater than 40 cm are intended for femoral vein insertion. This catheter is indicated for > 30 days long term placement. The valved peelable introducer is

intended for use in the percutaneous insertion of catheters in the venous system.

CONTRAINDICATIONS: The catheter is intended for long-term vascular access only and should not be used for any purpose other than indicated in these instructions. The valved peelable introducer sheath is not designed for use in the arterial system.

WARNINGS AND PRECAUTIONS: Please see package insert for complete list of warnings and precautions

POTENTIAL COMPLICATIONS: Consult package insert for a complete list of potential complications.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.



USA 14 Plaza Drive, Latham, NY 12110 tel: 800-772-6446 or 518-798-1215 fax: 518-798-1360
 International Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O The Netherlands
 tel: +31 (0)20 753 2949 fax: +31 (0)20 753 2939