


RESILIENCE™

Resilience Through-the-Scope Esophageal Stent




RESILIENCE™

Resilience Through-the-Scope Esophageal Stent



Anti-Migration Performance

Proximal and distal double flanges help to resist migration, anchoring the stent in place.^{1,2}



A Through-the-Scope Stent


Resilience decreases procedure time compared to the over-the-wire technique.³

Accurate, Easy, and Controlled Placement

Designed with a single-handed deployment system.¹


Safety Buttons

Prevent premature stent deployment and ensure proper alignment.¹



Unique Size Options¹

Only stent on the market offered in 14 mm x 50 mm, 17 mm x 50 mm, and 20 mm x 50 mm for short strictures.¹



14 mm x 50 mm

17 mm x 50 mm

20 mm x 50 mm

ORDERING INFORMATION

Catalog Number	Stent Midbody	Total Midbody Length	Flare
RE145	14 mm	50 mm/10 mm	24 mm
RE146	14 mm	60 mm/20 mm	24 mm
RE147	14 mm	70 mm/30 mm	24 mm
RE175	17 mm	50 mm/10 mm	27 mm
RE176	17 mm	60 mm/20 mm	27 mm
RE177	17 mm	70 mm/30 mm	27 mm
RE205	20 mm	50 mm/10 mm	28 mm
RE206	20 mm	60 mm/20 mm	28 mm
RE207	20 mm	70 mm/30 mm	28 mm

1. Data on File.
2. Mozafari et al. 2018. "Migration Resistance of Esophageal Stents: The Role of Stent Design." *Comput Biol Med* 1 (Sept): 43–49. doi: 10.1016/j.complbiomed.2018.06.031.
3. Mangiavillano et al. 2021. "Through-the-Scope Esophageal Stent for the Relief of Malignant Dysphagia: Results of a Multicentric Study (with video)." *Dig Endosc* 33, no. 1 (Jan): 118–24. doi: 10.1111/den.13704.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733