

C₂ CryoBalloon[®]

Ablation System Accessories

NEXT-LEVEL EFFICIENCY



C₂ CryoBalloo[®]

Ablation System Accessories



C₂ CryoBalloo NitroClip™



Interface device for expelling nitrous oxide from the C₂ CryoBalloo Controller into the hospital evacuation system.



C₂ CryoBalloo EndoGrip™



Optional device to reduce balloon migration. Once the catheter is in position, depress the EndoGrip and insert it onto the catheter, between the blue rubber grips, flush to the biopsy valve.



C₂ CryoBalloo SIDECAR™



External working channel if a therapeutic endoscope is not available. For use with a non-therapeutic endoscope.



C₂ CryoBalloo Storage Case



Hard-shell polypropylene exterior to protect equipment from inadvertent damage.

Holds two C₂ CryoBalloo Controllers and one C₂ CryoBalloo Foot Pedal.

ORDERING INFORMATION

Catalog Number	Description
C₂ CryoBalloo Accessories	
FG-1040	EndoGrip (5 Pack)
FG-1013	Cartridge (5 Pack)
FG-1054	Storage Case
FG-1011	SIDECAR
FG-1050	NitroClip
C₂ CryoBalloo Miscellaneous	
PRT-2434	Replacement Cap



C₂ CryoBalloo Cartridge (5 Pack)



C₂ CryoBalloo Replacement Cap

Attaches to the nitrous-oxide cylinder and secures to the C₂ CryoBalloo Controller.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733